

PEMS Training Manual

Date Created: April 30, 2004

Last Updated: July 27, 2011



ABSTRACT

This document will provide a step-by-step series of instructions that will guide a user in the creation, definition, and maintenance of a project in the Project Environmental Measurement System (PEMS). It will provide general information on the flow of data through PEMS and is not intended to be all-inclusive regarding day-to-day operations.

ACRONYMS

CAS	Chemical Abstracts Service
EDD	Electronic data deliverable
FCOC	Field chain of custody
FCTM	Field Collection Task Map
FPMC	Flow Paced Monthly Composite
LCOC	Laboratory chain of custody
OREIS	Oak Ridge Environmental Information System
PEMS	Project Environmental Management System
PC	Personal computer
QC	Quality control
RTL	Ready to load
SAP	Sample Analysis Plan
SOW	Statement of Work

1.0 Introduction and Access	6
2.0 Global Administration	8
3.0 Project Administration	10
3.1 Project Definition	12
3.2 Project Equipment	13
3.3 Project Lab	15
3.4 Project Parameter Analysis Group	16
3.5 Project Statement of Work	19
3.6 Project Task	23
3.7 Project Task Location	28
3.8 Project User Role	32
4.0 Data Management	33
4.1 Sample Analysis Plan	34
4.1.1 Transfer and Validate SAP Data	35
4.1.2 Get a Hardcopy of SAP Data from the Database	41
4.1.3 SAP Data File History	43
4.1.4 SAP Parameter Group Descriptions Report	44
4.1.5 SAP Parameter Group Definitions Report	46
4.1.6 SAP Sample Locations Report	48
4.1.7 SAP Comparison Report	50
4.2 Planning	51
4.2.1 Field Collection Task Reference	52
4.2.2 Field Collection Task Map	57
4.2.3 Field Chain of Custody (FCOC)	65
4.2.4 Flow Paced Monthly Composite (FMPC)	79
5.0 Field Data	88
5.1 Enter Collection Date and Time	88
5.2 Laboratory Chain of Custody (LCOC)	90
5.2.1 Printing LCOCs	94
5.2.2 Electronic Export LCOC	95
5.3 Bar Code Operations	97
5.3.1 Print Bar Code Menu Sheets	100
5.3.2 Transfer Barcode Reader Data	103
5.3.3 Edit Sample Barcode Reader Data	106
5.4 Field Measurements	112
5.5 Associate Samples	119
5.6 Associate Samples to Field Measurements	124
6.0 Laboratory Data	128
6.1 Transfer and Validate EDD Data	129
6.2 Get a Hardcopy of EDD Data from the Database	134
6.3 EDD Data File History	136
6.4 SOW Detail Information	138
6.5 Missing Analytes Report	139
7.0 Analytical Data	142
7.1 Validation Qualifier Assignment	142
7.2 Edit Sample Results	150
7.3 Apply Data Assessment Codes	151
7.4 Result/DRG Verification	155
8.0 RTL Review and Generation	159
8.1 RTL Review	159
8.2 Send a RTL to OREIS	161
I. CONTACTS	166
II. PROJECT INFORMATION	166
III. DATA INFORMATION (attach pages as needed)	166
IV. OREIS Receiving/Review	166

9.0 Utility Functions 168

 9.1 Delete/Rename Sample..... 168

 9.2 Change Sample Status 172

 9.3 Change Monitoring Location..... 173

10.0 Reports 174

1.0 Introduction and Access

Opening an Account

To open an account on PEMS, please send an e-mail to pems@ettp.doe.gov.

After obtaining a PEMS account, the user should have access to at least one project. The user will be assigned a user id and a password. Users' passwords will change every 6 months. Support personnel will notify all users of their new passwords. Figure 1 displays the log in screen.




Figure 1

Users may access projects by project name or by SOW number as shown in Figure 2.




Figure 2

2.0 Global Administration

Information for all Users

Figure 3 shows the various lists that may be accessed by all users. These lists contain the valid values currently in PEMS and OREIS used for all projects. Additions may be made only to the 'Update' option. Applications that are assigned the 'Read Only' level will only display information for the specified field. A badged person, monitoring location, preservative, sampling device, user information, company, parameter analysis group, sample container type, and sampling procedure may be added. Other additions or changes must be made by a change request to the PEMS, OREIS, or Tracker support personnel.




Figure 3

Some fields used in the PEMS system are shared by two other Environmental Information (EIM) systems. The applications listed under the ‘OREIS related tables’ header display information that is present in both the PEMS and OREIS systems. PEMS synchronizes the information for these fields twice daily to ensure that the latest information added to OREIS is included in PEMS. The applications listed under the ‘Tracker Information’ header list data that is present in the Tracker system. As this information is retrieved from Tracker when needed, there is no refresh for these fields.

The names of the applications and a brief summary of the information they display are presented in the following table.

Application Name	Purpose
Analysis Method	List of valid methods
Analysis Type	List of analysis types and descriptions
Badged Person	List of persons by badge number, user id and company
Company	List of companies managing projects
Data Assessment	List of the valid data assessment codes and their definitions
Database Role	List of roles the DBA can assign personnel
Data Form Type	List of forms that can be used. Default is FAS-01, field COC
Equipment Type	List of equipment that can be tracked to a low level
Field Activity Method	List of equipment that can be tracked to a low level
Field Equipment	List of equipment by id, includes active and inactive
Laboratory	List of laboratories that have or are performing work
Matrix Type	List of matrices by code
Measurement Unit	List of results units used by the laboratories
Monitoring Location	List of locations that can be assigned, can be by project
Parameter	List of analytes that can be requested by the project
Parameter Analysis Group	List of analyses that can be grouped together by container
Preservative	List of preservatives for containers
Sample Container Type	List of containers by size and material
Sample Status	List of statuses a sample can be assigned
Sample Type	List that defines quality control/assurance samples
Sampling Device	List of sampling devices that can be assigned
Sampling Procedure	List of sampling procedures that can be assigned and used
SOW	Tracker or non-Tracker statement of work number
User Information	User can update own information
Valid Qualifier Reason	List of validation qualifier reasons that can be assigned
Valid Result Qualifier	List of validation qualifiers that can be assigned to data

3.0 Project Administration

To be able to create a project, the user must have either the Project Data Coordinator role or the Sample Planner role. After logging in to the project, the user should select the Project Admin menu from the PEMS Main Menu, as indicated below in Figure 4.




Figure 4

Once Project Admin has been selected from the PEMS Main Menu, the Project Administration menu is displayed as shown in Figure 5. This menu contains all of the applications used for defining a project. From this menu, the project's laboratories, sampling locations, statements of work, user access and other information can be entered. The menu is arranged in the order of the information that should be entered. For example, a project's sampling locations by task cannot be entered until the project's task code is entered.




Figure 5

3.1 Project Definition

To create a new project, the user should select the Project Definition application from the Project Administration menu. This application will display the information for the current project as shown in Figure 6. Project information includes the project code, project name, project description, project manager and other data about the project.

At the top of the screen is a button labeled 'Add A Project'. Clicking on this button will begin the project creation process. The Project Definition screen will be re-displayed with the fields cleared for entry of a new project. Required fields are marked with arrows.

When all of the fields have been entered, click on the "Click Here To Continue" icon at the bottom of the screen. If the values are valid, a new project will be created, the user will be automatically switched to the new project and the next step of the project creation process will be displayed.

The screenshot shows a Netscape browser window titled "PEMS - Netscape". The address bar displays "https://www.pems.bechteljacobs.org/qpems/qpems/Pems.load_menus". The page features a blue sidebar menu on the left with the following items: "Christy's Test Project", "Lab Data", "Data", "Management", "Global Admin", "Project Admin" (highlighted in yellow), "Reports", "What's New", "Switch Project", and "Help". The main content area is titled "Project Definition" in large blue font, with a subtitle "Required fields are marked with" followed by a red arrow icon. Below the title is a button labeled "Add A Project". The form contains the following fields:

ID No.:	CHRISTY
Sponsor:	Whoever
Name:	Christy's Test Project
Begin Date(mm/dd/yyyy):	09/12/1998
Sample Org Name:	Christy's Sampling Team
OU Acronym:	dddd
RQMT Source ID:	dddd
End Date(mm/dd/yyyy):	12/12/1998
ESO Number:	
Manager Name:	
Manager Phone:	
Technical Contact:	
Technical Contact Phone:	
WBS No.:	
Description:	

At the bottom of the browser window, the status bar shows "Document: Done".

Figure 6

3.2 Project Equipment

Once the new project has been created, the user will be presented with the Project Equipment screen, as shown in Figure 7. This application is used to provide tracking for equipment that will be used by the project.




Figure 7

To assign equipment to the project, the user should click on the 'Change Field Equipment For This Project' button at the top of the screen. This button will display a list of all equipment that is currently available for the project to use (Figure 8). This list of equipment is maintained as part of the PEMS global support tables.

If the equipment to be used is not listed, the user must go to the Global Administration menu (from the PEMS Main Menu) and select the Project Equipment application to add the new equipment. Once the

equipment has been added to the appropriate support table, the user should return to the Project Equipment application to assign the equipment to the project.




Figure 8

To add equipment to the project, the user should click on the check box beside each individual piece of equipment. When all of the equipment to be added has been selected, click on the 'Click Here To Continue' icon to return to the Project Equipment main screen (Figure 7). Equipment can be added to a project at any time and having project equipment defined is not mandatory for project creation.

3.3 Project Lab

To continue with the project creation process, click on the 'Click Here To Continue' icon from the Project Equipment main screen. This will display the Project Lab main screen (Figure 9). In the past, this application was used to assign laboratories to the project. Assignment of laboratories is now performed by the Project Statement of Work application. This application has been retained to give a quick way to display which laboratories are assigned to the project.




Figure 9

3.4 Project Parameter Analysis Group

To continue with the project creation process, the user should click on the ‘Click Here To Continue’ icon, which will display the Project Parameter Analysis Group main screen. This screen is shown in Figure 10.




Figure 10

It is likely that a code representing the analyses that will be performed has already been created. To display the current list of analytical codes, the user should click on the ‘Change Analysis Group(s) For This Project’ button. Depending upon the configuration of the user’s PC and/or the user’s network, this screen may take additional time to load. PEMS will then display a new screen (see Figure 11) that shows all of the Analytical Group codes defined in the system.

To assign the analytical group code(s) to the project, click the check box next to the desired code and enter a value in the 'Label Qty' field for that code. The 'Label Qty' column is used for the printing of sample container labels. Entering a value of '3', for example, in this column would cause 3 sample container labels to be printed for any sample that is assigned to this analytical group. Depending upon the resolution of the user's PC, this column may be off-screen. If this is the case, a scroll bar will appear at the bottom of the browser window. Use this bar to display this field so that the label quantity can be specified.

Alternately, the user may use the Sample Analysis Plan (SAP) load to enter information. (See 4.1) If the user has assigned some analytical group codes to the project and no longer wishes them to use them, the analytical group codes can be removed from the project using this screen. However, it is only possible to remove an analysis group code if that code has not been used in any sample plan. The analytical group code can be unassigned from the project by de-selecting the check box and removing the quantity value from this screen. Clicking the 'Click Here to Continue' icon will save the changes.

PEMS - Netscape

File Edit View Go Communicator Help

Back Forward Reload Home Search Netscape Print Security Shop Stop

Bookmarks Location: https://www-pems.bechteljacobs.org/qpems/qpems/Pems.load_menus What's Related

Instant Message Internet EIM Sites Lookup New&Cool LMES Internal H Oracle Document CNET.com - Welc WebWHOS: Search Seagull's Therm

Add Parameter Analysis Group To Project

CLICK HERE TO CONTINUE

	Code	Name	Label Qty
<input type="checkbox"/>	%SOLIDS	Percent Solids	<input type="text"/>
<input type="checkbox"/>	(BTU)	BTU	<input type="text"/>
<input type="checkbox"/>	0123456789	0123456789012345678901234567890123456789012345678901234567890123456789012345678901234567890123456789	<input type="text"/>
<input type="checkbox"/>	1194RAD1	Radioactive Set 1 (BXPT1194)	<input type="text"/>
<input type="checkbox"/>	1194RAD2	Radioactive Set 2 (BXPT01194)	<input type="text"/>
<input type="checkbox"/>	1227CHEM1	SVOC, Pesticides, Ignitability, pH, Oil & Grease, Paint Filter	<input type="text"/>
<input type="checkbox"/>	1227CHEM2	SVOC, Pesticides, Herbicides (Water)	<input type="text"/>
<input type="checkbox"/>	1377RADS	All SOW RADs	<input type="text"/>
<input type="checkbox"/>	1501-1	NIOSH 1501, BENZENE ONLY	<input type="text"/>
<input type="checkbox"/>	1501-2	NIOSH 1501- Toluene, Xylene, Benzene	<input type="text"/>
<input type="checkbox"/>	1501-3	NIOSH 1501 - Toluene & Xylene	<input type="text"/>
<input type="checkbox"/>	1514RADS	Am-Iso,Np-237,Pu-Iso,Th-Iso,U-Iso,%U-235>Total U,TC-99	<input type="text"/>
<input type="checkbox"/>	1566RADS	Am-241&Np-237&Pu-238,239/240&Tc-99&U-233/234,235,236,238&%U-235&TotalU	<input type="text"/>
<input type="checkbox"/>	1570RADS	Am-Iso,Pu-Iso,Th-Iso,U-Iso,Tc-99,Sr-90,Np-237,TotalU,Gamma Spectroscopy/PER SOW	<input type="text"/>

Document: Done

Figure 11

When the edit process is finished, click on the 'Click Here To Continue' icon to save the selections and assign the analytical groups to the current project. The display will return the current listing of analytical groups for the project. Additional analysis group codes can be added to the project later, if needed. To continue with project definition, click on the 'Click Here To Continue' icon, which will display the Project Statement of Work (SOW) application (see Figure 12).

3.5 Project Statement of Work

The Project SOW application creates statements of work that are issued to labs for the analysis performed on collected samples and edits information regarding these statements of work. The main screen of this application shows a drop-down list that currently has 2 options, 'Add New SOW From Tracker' and 'Add A SOW That Is Not In Tracker'. An existing SOW could be edited by selecting the SOW from the drop-down list and clicking on the 'View Details' button.

If the Tracker option is selected, then only SOWs that are defined in the Tracker system will be available for creation in PEMS. Selecting the SOW will also define the laboratory for that SOW and assign that laboratory to the project.




Figure 12

If the Tracker option is selected, a prompt will be displayed in a pop-up window (Figure 13) asking for the numeric SOW to add to PEMS. Enter a numeric value and click the 'OK' button.




Figure 13

If there are multiple SOWs that match the numeric value entered, the pop-up window will be replaced with a second pop-up window (Figure 14) showing all of the statements of work that match that value in a drop-down list. Select the appropriate SOW and click on the 'OK' button to create the statement of work. If no statements of work match the entered numeric value, the drop-down list will be empty and the 'OK' button will not be present on the pop-up window. Click on the 'Back' button to try a different numeric value.




Figure 14

If the entered numeric value matches only one SOW in the Tracker system, the SOW will be created and the lab will be assigned to the project. If the lab associated with the entered SOW does not exist in PEMS, the lab code and lab name will be retrieved from Tracker and inserted into PEMS. An edit screen will then be displayed so that the user can enter additional information about the SOW. This

screen is the same screen that will be displayed if the user edits an existing SOW. On creation of a new Tracker SOW, the appearance of this edit screen may vary depending on certain conditions. If the lab associated with the Tracker SOW does not have a PEMS login account, an informational message will also be displayed along with the edit screen, as shown in Figure 15.

This message will instruct the user that, though the SOW has been created, the lab will not be able to load analytical results because they do not have a PEMS account. It will give instruction to the user on how the lab can acquire an account by linking to information on the PEMS Home Page. The user should direct the laboratory to request a PEMS account following the procedures outlined in the "Getting Started" section of the PEMS Home Page. This message is only displayed at creation time of the SOW. It will not be displayed if the SOW is edited later.

The screenshot shows a Netscape browser window with the URL https://www-pems.bechteljacobs.org/qpems/qpems/Pems.load_menus. The page title is "SOW Definition". On the left is a blue navigation menu with links: Instruction, Manual Project, Lab Data, Data, Management, Global Admin, Project Admin (highlighted), Reports, What's New, Switch Project, and Help. The main content area contains a form with the following fields:

ID No.	INSTMAN
SOW No.	BKET01320
SOW Rev No.	0
Effective Date(mm/dd/yyyy)	10/07/1999
Lab Code	DCSL (DCSL)
Effective End Date(mm/dd/yyyy)	<input type="text"/>
BJC Tracker SOW?	Y
Fund Amount	<input type="text"/>
Required No.	<input type="text"/>
Charge No.	KKUPH4F3
RQMT Source ID	<input type="text"/>
Tracking No.	<input type="text"/>

Below the form is a red button labeled "CLICK HERE TO CONTINUE".

Figure 15

None of the fields on the additional information screen are required. This information may be added or changed at the user's discretion. When finished, clicking on the 'Click Here to Continue' icon will return the user to the main screen of the SOW application where more SOWs may be created, if needed.

If Tracker will not be used to define the SOWs, the user will be able to specify any laboratory defined in the PEMS system and any SOW number. In this case, the format of the SOW number need not follow the Tracker specification for SOWs.

Upon selecting this option, a screen will be displayed that prompts the user for the SOW number and shows a drop-down list containing all of the labs currently defined in PEMS (See Figure 16). The lab must be present in the PEMS system to be selected.

The screenshot shows a Netscape browser window titled "PEMS - Netscape". The address bar displays the URL "https://www.pems.bechteljacobs.org/qpems/qpems/Pems.load_menus". The left sidebar contains a navigation menu with the following items: "Instruction Manual Project", "Lab Data", "Data Management", "Global Admin", "Project Admin" (highlighted in yellow), "Reports", "What's New", "Switch Project", and "Help". The main content area has the heading "Add A SOW That Is Not In Tracker" in large blue font. Below the heading, it states "Required fields are marked with" followed by a red circle with a plus sign. A sub-note says "If the desired laboratory is not present in this list, Contact PEMS Support". The form consists of two fields: "SOW No.:" with a text input box and a red circle with a plus sign, and "Laboratory:" with a dropdown menu showing "Babcock and Wilcox, Lynch-BWLVA". Below the form is a red button labeled "CLICK HERE TO CONTINUE". The status bar at the bottom shows "Document: Done".

Figure 16

After entering the SOW number and selecting the lab, the same additional edit screen will be displayed as in the Tracker SOW. As in the Tracker SOW, if that lab does not have a login account for the PEMS system, a message will be displayed providing information on acquiring a PEMS account. When finished with data entry, clicking on the 'Click Here to Continue' icon will return the user to the main screen of the SOW application.

3.6 Project Task

From the main screen of the SOW application, clicking on the 'Click Here to Continue' icon will continue the project definition process by displaying the Project Task application. The main screen of this application is shown in Figure 17. Since this is a new project, there are no tasks defined.




Figure 17

A project task is a 5-character user-defined code that describes a level of activity in the project. This activity can be any sampling activity that the user wishes to model in PEMS. For example, it could be the sampling activity for all groundwater wells for the entire year. It could also represent sampling activity for one well for one occurrence. The definition of the task is completely at the user's discretion.

Two buttons will be displayed at the top of this screen. The first one, 'Add New Task(s) To This Project', will display a task definition screen as shown in Figure 18.

PEMS - Netscape

File Edit View Go Communicator Help

Back Forward Reload Home Search Netscape Print Security Shop Stop

Bookmarks Location: https://www.pems.bechteljacobs.org/qpems/qpems/Pems.load_menus What's Related

Instant Message Internet EIM Sites Lookup New&Cool LMES Internal H Oracle Document CNET.com - Welc WebWHOS: Search Seagull's Therm

Instruction Manual Project

Lab Data

Data

Management

Global Admin

Project Admin

Reports

What's New

Switch Project

[Help](#)

Add Task To Project

Copy The Tasks From An Existing Project

Code	Description
<input type="text"/>	<input type="text"/>
<input type="text"/>	<input type="text"/>
<input type="text"/>	<input type="text"/>
<input type="text"/>	<input type="text"/>
<input type="text"/>	<input type="text"/>
<input type="text"/>	<input type="text"/>
<input type="text"/>	<input type="text"/>
<input type="text"/>	<input type="text"/>
<input type="text"/>	<input type="text"/>
<input type="text"/>	<input type="text"/>

CLICK HERE TO CONTINUE

Document: Done

Figure 18

The task definition screen allows for the creation of new tasks for the project. The user should enter the task code and a description of the activity represented by this task. The entered task code must be unique for the project. When data entry is completed, click on the 'Click Here to Continue' icon to return to the main screen of this application.

A button, labeled 'Copy The Tasks From An Existing Project' is displayed at the top of the Add Project Task screen. This button permits data entry by allowing the user to copy the tasks from another project. Some users create new projects every year to reflect that year's activity, but will be performing the same tasks in that new project as used in a previous project. Clicking on the button will display a list of all of the projects that the user has privileges to access (See Figure 19). This screen is very similar to the 'Login by Project' screen that is displayed during the login process.

ID	Name
BIOASSAY	BJC Radworker Bioassay Sampling
BJC09914	PORTSMOUTH - SW-3 FLOOR SWEEPINGS
BJC09914A	PORTSMOUTH - SW-3 FLOOR SWEEPINGS
CHRISTY	Christy's Test Project
CNF 2002	2002 CNF NPDES Monitoring
CNF INFLUE	Radian Influent Sampling
DD356	Donna's Test Project
DD357	Donna's second test Project
EH INVESTI	EH Investigation ETPP May 2000
ETTP-MISC	ETTP Miscellaneous Sampling Events by MDM
FLOW-PACED	Test of Task Ref Tbl. and Flow Paced Composites
IWQP-COMP	Integrated Water Quality Program - Composite
IWQP-ETTP	ETTP Integrated Water Quality Program (FY 1999)
K201	Stormwater Pollution Prevention Program
NANCY	Nancy Sampling at 1035
PGDP-99	Portsmouth Gaseous Diffusion Plant 1999
PORTSMAY02	PEMS Training
QA990118	qa project
RW8	Semi Annual
SCOTTSPRJ	training test

Figure 19

Clicking on the project code will display an intermediate screen showing all of the tasks defined for the selected project (See Figure 20). A check box is displayed beside each task code. To copy the task to the new project, click on the check box to select the code. The user may select any or all of the task codes to transfer them to the new project. Clicking on the 'OK' button will copy the selected tasks and return the display to the main screen of the Project Task application.




Figure 20

The second button on Figure 21, 'Delete Task(s) From This Project', allows the user to remove tasks from the project. To delete a task, click on the check box next to the task code and definition to select it. Then, click on the 'Click Here to Continue' icon to remove the task from the project. The display will return to the main screen of the application. Please note that the 'Delete Task' function will not allow the deletion of a task that is currently in use by the project on existing documentation (FCOC, LCOC).




Figure 21

3.7 Project Task Location

Clicking on the 'Click Here to Continue' icon from the Project Task main screen will display the final step of the project definition process. This step is the Project Task Location application, which is displayed in Figure 22. Since this is a new project, there are no locations assigned to the tasks created in the previous step.




Figure 22

To assign locations to the tasks defined, click on the 'Change Task Locations For This Project' button. This will display an intermediate screen that shows two drop-down lists. The first list contains all of the tasks defined for the project. The second list is a group selector for all of the locations defined in the PEMS system. This second list is used to limit the locations displayed to those assigned to a specific site to speed data entry. For multi-site projects, an option is present to show all locations. This option is the default.




Figure 23

To assign locations to the tasks, select the appropriate task and site and click on the 'Click Here to Continue' icon. This will select the task and display a list of all locations available for the select site option (See Figure 24).




Figure 24

A check box is present beside each of the listed locations and their descriptions. To add a location to a task, click on the check box to select the location. If the intended location is not present in the list, it may need to be added to the PEMS system through the Monitoring Location application under the Global Administration menu. The user should confirm that this location is not in the PEMS system before adding a new location. Contact PEMS Support if assistance is needed in finding a location. When finished with the selection of locations, click on the 'Click Here to Continue' icon to add the locations to the task and return to the main screen of the application.

The same locations may be used for different tasks, as often as needed. If editing an existing task location list, simply un-check the checkbox of a task location to remove that location from the task. However, once a location and task code combination have been used in a field chain of custody, that location cannot be deleted from the task. No message will be displayed if the user attempts to delete the location. The operation simply will not succeed.

The definition of project is now complete. The user may now begin planning the sampling efforts. How the user proceeds may depend upon the needs of the project. At this point, a user may go straight to the Field Chain of Custody application to create sample planning documentation. If the project will need more regular sampling or has an extensive volume of samples, it is advised that the user set up a formal sample plan using the Task Reference and Field Task Map applications.

3.8 Project User Role

PEMS has a variety of user 'roles' which allow project personnel different levels of access to the project information. (Figure 25) Clicking on the appropriate box will give each user that role of access for each project.

The roles are:

- A. Data Validator/Verifier – Allows assignment of validation or verification qualifiers
- B. Lab – Allows labs to load data
- C. Project Data Coordinator – Allows full project access
- D. Reports Viewer – Allows 'read only' access to the PEMS Report and Global menus
- E. Sample Planner – Allows access to sample planning functions




Figure 25

4.0 Data Management

To begin the creation of a sample plan, the user should first click on the 'Data Management' menu option on the left pane from the PEMS Main Menu. This will display the Data Management menu, as shown in Figure 26. This menu is where the majority of the day-to-day PEMS applications can be found.




Figure 26

4.1 Sample Analysis Plan

Comma-separated value (CSV) files taken from a Sample Analysis Plan (SAP) (Figure 27) may be inserted into tables in PEMS when creating a project in PEMS. The user must have either Project Data Coordinator (PDC) or Sample Planner access to the current project to be able to load the SAP files.




Figure 27

4.1.1 Transfer and Validate SAP Data

The file structures for each of the SAP files are listed below. The column header row should be removed from the file before loading into PEMS.

File 1: Parameter Group Descriptions (comma delimited). The file should be saved in text format and file the name should end with .des.

Field Number	Field Name	Max Width	Definition
1	FY	4	Fiscal Year
2	Parameter_Group	10	Parameter Group code
3	Parameter_Group_Description	80	Parameter Group description from PEMS

The second file is Parameter Group Definitions (comma delimited). This text file should end in .def.

Field Number	Field Name	Max Width	Definition	Valid Values
1	FY	4	Fiscal Year	
2	Parameter_Group	10	Parameter Group code	See File 1 structure: Parameter_Group field for list of valid Values
3	Analyte_Name	90	Chemical name for the analyte	See Chemical Name in OREIS Parameter Code list
4	Detection_Limit	14	Method Detection Limit	
5	Detection_Unit	10	See Units in OREIS Unit list	
6	Analyte ID	9	CAS number or pseudo-CAS number. This field is optional	See Parameter in OREIS Parameter Code list

The third file defines sample locations. This file should end with .smp.

Field Number	Field Name	Max Width	Definition	Valid Values
1	FY	4	Fiscal year	
2	Project_ID	10	Identifier for the project the sample is a part of	Must match an existing PEMS project ID
3	Watershed_Name	10	Watershed location	
4	Task	5	Task Grouping	
5	Station_Name	15	Sample Location	See PEMS Monitoring Locations in Global Admin Section
6	Matrix_ID	2	A code for the sample matrix	See OREIS MED-TYPE
7	Sampling_Frequency	20	Sampling Frequency	WK – weekly, M – Monthly, Qx fiscal quarter
8	Field_Duplicate	1	Is Field Duplicate required	Y or N
9	Parameter_Group	100	List of parameter groups	See File 1 structure: Parameter_Group field for list of valid values

These files may be loaded into the PEMS database using the 'SAP Data Load' function. The new project must exist in PEMS before the .des and .def files can be loaded into PEMS. SAP files must be unique for the project.




Figure 28

When each file is loaded, the system will run a set of checks on the data. (Figures 29-30)




Figure 29

The system will confirm that the data has been inserted into the database.




Figure 30

The system will assign a unique file identifier to each file as shown in Figure 31. All SAP files will have an identifier that starts with "S." The SAP files are loaded into the current project. Creation of a project (Project Menu: Project Definition) must have already occurred to receive the files that define the sample plan.




Figure 31

The user may now load the files into the new project. The .smp file loads the Field Collection Task Reference screen in the new project. The Tracker user would then pull these files from PEMS to create portions of the Tracker SOW.

4.1.2 Get a Hardcopy of SAP Data from the Database

Users may obtain a copy of the SAP data by date range, for the current or any project to which the user has access by using the 'Get A Hardcopy of SAP Data From The Database' function.




Figure 32

Figure 33 shows an example of the SAP file information displayed by this function.




Figure 33

4.1.3 SAP Data File History

The user may enter the database number for a particular SAP file and view information about the file (e.g., when it was loaded). Files can be downloaded through the Hardcopy option and viewed locally through a local editor. The file prompt will accept a file identifier 'Sxxxxx' or an actual name (e.g., sapdes.def).




Figure 34

4.1.4 SAP Parameter Group Descriptions Report

The 'SAP Parameter Group Descriptions Report' gives a list of all the group codes with corresponding descriptions found in PEMS for the selected fiscal year. (Figure 35)





Figure 35

Figure 36 shows the SAP Parameter Group Descriptions Report in a printable .rtf format.



Fiscal Year	Parameter Group Code	Parameter Group Description
2004	SO-GAS-04	Soil Gas VOCs
2004	SO-HG-04	mercury
2004	SO-RAD1-04	U-Isotopic;Gross Alpha/Beta
2004	SO-RAD2-04	Sr-90;U-Isotopic;Gross Alpha/Beta;Gamma Spec
2004	SO-RAD3-04	Sr-90;U-Isotopic;Gross Alpha/Beta;Gamma Spec;Cm-243/244
2004	SO-VOC-04	VOCs
2004	W-BIO1-04	Alk.;BOD;Nitrate;Nitrite;Cl;F;Sulf
2004	W-BIO2-04	Ammonia;COD
2004	W-BIO3-04	Sulfide
2004	W-BIO4-04	Methane;Ethane;Ethene
2004	W-BIO5-04	TOC
2004	W-BIO6-04	Fe_2
2004	W-BIO7-04	Fe_3;Mn
2004	W-HG-04	mercury
2004	W-MET-04	Dissolved Metals
2004	W-RAD1-04	Gross Alpha/Beta;Gamma Spec
2004	W-RAD2-04	U-Isotopic;Gross Alpha/Beta
2004	W-RAD3-04	Sr-90;U-Isotopic;Gross Alpha/Beta;Gamma Spec

Figure 36

4.1.5 SAP Parameter Group Definitions Report

The SAP Parameter Group Definitions Report will give a list of analytes, limits, units and Chemical Abstract Numbers (CAS) associated with each group for the year.




Figure 37

Figure 38 shows the SAP Group Definitions Report in a printable, .rtf format.

VUVA4GJ9.RTF - Microsoft Word

File Edit View Insert Format Tools Table Window Help

100% Times New Roman 9 B I U

BECHTEL JACOBS
Bechtel Jacobs Company LLC

Project Environmental Measurement System (PEMS)
SAP Parameter Group Definitions

Page: 1
Date: 04/21/2004 8:32 AM

Fiscal Year	Parameter Group Code	Analyte Name	Analyte ID	Detection Limit	Detection Unit
2004	SO-GAS-04	1,1-Dichloroethene	75354	25.0000	NG
2004	SO-GAS-04	1,2-Dichloroethene	540590	25.0000	NG
2004	SO-GAS-04	Benzene	71432	25.0000	NG
2004	SO-GAS-04	Tetrachloroethene	127184	25.0000	NG
2004	SO-GAS-04	Toluene	108883	25.0000	NG
2004	SO-GAS-04	Trichloroethene	79016	25.0000	NG
2004	SO-GAS-04	Vinyl chloride	73014	50.0000	NG
2004	SO-HG-04	Mercury	7439976	0.0500	MG/KG
2004	SO-RAD1-04	Alpha activity	12587461	5.0000	PCI/G
2004	SO-RAD1-04	Beta activity	12587472	5.0000	PCI/G

Draw AutoShapes

Page 1 Sec 1 1/7 At 1.9" Ln 1 Col 1 REC TRK EXT OVR

Start 8:33 AM

Figure 38

4.1.6 SAP Sample Locations Report

The 'SAP Location Report' gives location, station, task, matrix and other information by fiscal year.




Figure 39

Figure 40 shows the SAP Sample Location Report for the fiscal year in a printable .rtf format.




Figure 40

4.1.7 SAP Comparison Report

A Project Data Coordinator or Sample Planner uses the 'SAP Comparison Report' after the laboratories have loaded electronic data deliverables into PEMS. This report compares what was asked for in the SAP (number of samples requested by analysis group, station, matrix, etc.) with results received. (Figure 41)




Figure 41

4.2 Planning

There are four subsections to the Planning menu. These sections are titled 'Field Collection Task Reference', 'Field Collection Task Map', 'Field Chain of Custody' and 'Flow-Paced Monthly Composite'. The Automatic Storm Event is currently unavailable.




Figure 42

4.2.1 Field Collection Task Reference

Click on the Field Collection Task Reference option to display the main screen of the Task Reference application, as shown in Figure 43. Depending upon resolution of the user's PC, portions of this screen may not be visible. If this is the case, a scroll bar is displayed at the bottom of the screen to allow the user to see all of the data.




Figure 43

The Task Reference application provides a means for the user to create a template for sample planning activity. This template can be used repeatedly to define samples for collection.

To begin, click on the 'New Record' button at the top of the screen. This will display a data entry screen for the Task Reference application (See Figure 44). The first field, Field Sequence, is a numeric field

that is intended to assist the user in ordering planned samples. For example, if the user wanted to group all Trip Blank samples on the task reference screen, a sequence number in the 500 range could be used for each entry (e.g. 500, 501, 502, etc.). All Regular samples could be ordered into the 1000 series and so on. The sequence number on the main screen orders the Task Reference entries. The Field Sequence number is mandatory and must be unique for each entry.

The screenshot shows a Netscape browser window displaying the 'Field Collection Task Reference' form. The form is titled 'Field Collection Task Reference' and includes a note: 'Required fields are marked with *'. There are two buttons at the top: 'Save Data' and 'Save Data and Add Another Task'. The form fields are as follows:

Field Sequence:	<input type="text"/>
Project Task Code/Location:	1407J/26-005
Multiple Locations?:	No
Matrix:	AA-Ambient Air
Collection Device:	20-2020 Photoionization air monitor-PID
Collection FOP/Revision:	—/Rev. 0 (No sampling procedure)
Container Volume:	<input type="text"/>
Container Unit:	<input type="text"/> Search
Container Type:	<input type="text"/>
Collection Method:	<input type="text"/>
Averaging Period:	<input type="text"/>
Sample Frequency:	1
Preservative:	— (NO PRESERVATIVE REQUIRED)
Analytical Group:	1227CHEM1 (SVOC, Pesticides, Ignitability, pH, Oil & Grease, Paint Filter)
Sample Type:	FB (FIELD BLANK)

Figure 44

Using the Project Task Code/Locations drop-down list, select the appropriate task and location combination. There is a special field labeled 'Multiple Locations?' that can have a value of 'Yes' or 'No'. When this field has a value of 'No' the sample defined will be for the task and location specified in the previous field. When it has a value of 'Yes', the same sample definition will be used for ALL locations associated with that task. This is a quick way of defining samples for a task if the same type of sample will be collected at multiple locations. Task Reference entries that are so defined will have a value of 'MULTI' for the Monitoring Location field on the main screen of the Task Reference application.

Clicking on the 'Search' button will open a pop-up window that displays all of the possible containers in the PEMS system as seen in Figure 45. These values are displayed in a drop-down list. Selecting the

container volume, unit and type and clicking on the 'OK' button will close the window and populate these fields on the Task Reference record.

For detailed descriptions of the remaining fields on this screen, please refer to the *PEMS Data Dictionary* and/or the *PEMS User's Manual*, on the PEMS Home Page, <https://www-pems.ettp.energy.gov/pems/home/home.html>.

When data entry is completed, click on the 'Save Data' button to save the record and return to the main screen of the Task Reference application. To remain on the edit screen and add additional Task Reference records, click on the 'Save Data and Add Another Task' button. This will save the current information and clear the Field Sequence field. The entries from the previous record will remain, as this helps to speed data entry for similar samples.




Figure 45

On return to the main screen of the Task Reference application, all of the entered task reference records are displayed on screen. To edit a task reference record, simply click on any field of that record and the record editor will be invoked for that record. The value in the Field Sequence field cannot be changed and there is an additional button at the top right of the screen. This button, labeled 'Delete This Task', will remove the record from the Task Reference template. It will not affect the tasks defined for the project or any of the project's samples already in the PEMS system.

After the task reference records have been created, users may speed data entry by using the duplicate mode to copy information. Clicking on the 'Duplicate Record' button will place the application in duplicate mode, which is indicated by an arrow icon next to the button (see Figure 46). To copy a record, click on any field of the record to be duplicated. To ensure that the entered sequence number is unique, make the desired changes to the new record and click on the save button to save the changes. The display will return to the Task Reference main screen.




Figure 46

After creating the Task Reference template records, the user may wish to re-order the records so that they are visually displayed differently. Clicking on the 'Renumber Sequence' button will display an edit screen that contains all of the records in the template (See Figure 31). On this screen, all of the template fields are display-only, with the exception of the sequence number field. By changing the values of the sequence fields, the user may re-group or re-sort the records in any order.




Figure 47

As before, each field sequence number value must be unique. Click on the 'Save' button to save the template and return to the Task Reference main screen. The Task Reference template may be modified at any time during the course of the project to add or remove records as needed. The user may set up all of the tasks for a project before continuing or set up just one task and continue on to the next phase of sample planning.

4.2.2 Field Collection Task Map

The Field Collection Task Map (FCTM) works with the template data stored in the Task Reference to create a more detailed sample plan. To access the FCTM, click on the 'Data Management' menu option in the left pane of the PEMS application to display the Data Management submenu. Click on the 'Field Collection Task Map' option to display the main screen of the Field Task Map, as shown in Figure 48.




Figure 48

The main screen of the FCTM shows a drop-down list and 3 radio buttons. The drop-down list defaults to the 'Create A New FCTM' option and the radio button defaults to the edit option. The drop-down list contains all of the field task maps created for the project that have not been loaded. Since this is a new project, the list is empty. The radio buttons direct what action should be taken to the task map displayed in the drop-down list. The Delete FCTM radio button will delete the selected task map.

The Load FCTM radio button will analyze the selected task map to determine if it is valid and if so, will load the information in the task map into PEMS. The Edit FCTM radio button will invoke the field task map editor for the selected task map.

To create a task map, click on the 'Click Here to Continue' icon with the 'Create A New FCTM' option and the Edit FCTM radio button selected. Since this is a new task map, the editor will not be displayed. This screen displays a drop-down list, 2 check boxes and a field for data entry. The drop-down list contains all of the task codes and matrix combinations that are recorded in the Task Reference template. The data entry field is used to enter the planned date that the sampling for the selected task code/matrix combination will be performed. This value for this field defaults to the current date.

The screenshot shows a Netscape browser window titled "PEMS - Netscape" with the address bar displaying "https://www.pems.bechteljacobs.org/qpems/qpems/Pems.load_menus". The browser's menu bar includes File, Edit, View, Go, Communicator, and Help. The toolbar contains icons for Back, Forward, Reload, Home, Search, Netscape, Print, Security, Shop, and Stop. The bookmarks bar shows various links like Instant Message, Internet, EIM Sites, Lookup, New&Cool, LMES Internal H, Oracle Document, CNET.com, WebWHOS, Search, and Seagull's Therm. The main content area is titled "Field Collection Task Map" and includes a note: "Required fields are marked with *". The form itself has a green header with two columns: "Task/Matrix" and "Sample Date (mm/dd/yyyy)". Under "Task/Matrix", there is a dropdown menu showing "1407JAWD". Under "Sample Date", there is a text input field containing "01/28/2003" and a red asterisk icon. Below the header, there are two checkboxes: "Select all locations defined for this task" and "Use automatic Sample Id generation". At the bottom of the form, there is a red button labeled "CLICK HERE TO CONTINUE". A mouse cursor is pointing at the button. On the left side of the browser window, there is a vertical blue sidebar with a menu containing: Instruction, Manual Project, Lab Data, Data, Management (highlighted), Global Admin, Project Admin, Reports, What's New, Switch Project, and Help.

Figure 49

The first check box, 'Select all locations defined for this task', will automatically prepare a sample entry on the task map for every location that is defined for the task code. If this box is checked, the system will look at the locations defined for the task/matrix combination in the task reference template and use

that information to build the sampling locations. If this box is not checked, the system will look at the locations defined for the task and display that information on another screen, as seen on Figure 50. This screen will list the monitoring location id and description with a check box to the left of this data. By clicking on the check box(es) of the desired location(s) and then on the 'Click Here to Continue' icon, the field task map will be created.

The second check box, 'Use automatic Sample ID generation', is used to speed data entry by having the system generate the sample container identifiers (hereafter referred to as sample ids). If this box is checked, the system will generate sample ids for each sample record in the task map. These sample ids will be of the form 'SMPLxxxxxx', where xxxxxx is a unique number generated by the system. This option is typically used for high-volume sample plans where the value used for the sample id has no formatting specifications. In the initial display of the task map editor, the generated sample ids may not be changed. If the user would like to implement some combination of generated and non-generated sample ids, the sample ids may be changed by returning to the main menu of the Field Task Map application and re-editing the task map.




Figure 50

Once the information has been specified, the field task map is created and an edit screen is displayed as shown in Figure 51. None of the data entry fields are mandatory, as the task map may be created and edited at any time. The Sample ID field should contain the identifier of the container that will hold the collected sample. It is a 15-character, case-sensitive field that can be any alphanumeric character, as well as the following special characters: #, /, - and .(period). Every record on the task map must have a sample id before it can be transferred from the task map into the PEMS system. Since the task map is essentially a “what-if” document, the same sample id may be used on multiple task maps, but not inside the same task map. When task maps are loaded, however, a check is performed to determine if the new sample id(s) have already been used for the project. If so, the load is aborted and an error message is displayed to the user.

The Associated Sample ID field is mandatory only under certain conditions. For example, if a sample is designated as a Field Replicate (sample type of ‘FR’), it must be associated with a regular sample (sample type of ‘REG’). The sample id of the regular sample is entered into the Associated Sample ID field. A check is performed during the loading process to ensure that this field is populated when necessary.

The FCOC# field is used to record the form number of the Field Chain of Custody (FCOC) that will be used by the sample collectors in the field. It is a 15-character, case-sensitive field that will accept any alphanumeric character. This field can be used to separate the samples onto different physical documents. For example, specifying samples X and Y on FCOC# 1 and sample Z on FCOC# 2 will create 2 different documents for collecting these 3 samples.




Figure 51

The Form Type field is a drop-down list that contains the document types supported by the PEMS system. Currently, PEMS supports the FCOC form and the Flow Paced Monthly Composite (FPMC) form. If a sample's document type is set to 'FCOC', that sample will be transferred to a FCOC document when the task map is loaded. The document can only be edited through the Field Chain of Custody application. If the document type is set to 'Composite', the sample will be transferred to an FPMC form when the task map is loaded. The Flow Paced Monthly Composite application must be used to edit the FPMC documents.

Unwanted samples on the task map may be removed by clicking on the 'Delete Task' check box, then clicking on the 'Delete Task(s)' button. Additional samples may be added by clicking on the 'Add Task' button, which displays an edit screen as shown in Figure 52. Since the task and the matrix are used in the definition of the task map, these fields are displayed, but cannot be edited. The other fields are the same as used for the task reference template, with the exception that the project location field is limited to those locations defined for the task. When data entry is complete, click on the 'Save Task' button to return to the task map editor. Note that the Sequence Number field has been automatically

supplied by the system. Any new sample added to a task map will be given a sequence number that is equal to the current maximum sequence number on the task map + 1. This is done to ensure uniqueness for each entry in the task map.

The screenshot shows a Netscape browser window with the URL https://www-pems.bchteljacobs.org/qpens/qpens/Pems.load_menus. The page title is "Field Collection Task Map". On the left is a blue sidebar with navigation links: Instruction, Manual Project, Lab Data, Data, Management, Global Admin, Project Admin, Reports, What's New, Switch Project, and Help. The main content area contains a form with the following fields:

Project Task Code:	1407J
Project Location:	26-005
Matrix:	WG-Groundwater
Collection Device:	20-2020 Photoionization air monitor-PID
Collection FOP/Revision:	—/Rev.0 (No sampling procedure)
Container Volume:	
Container Unit:	<input type="text"/> Search
Container Type:	
Collection Method:	
Averaging Period:	
Preservative:	— (NO PRESERVATIVE REQUIRED)
Analytical Group:	1227CHEM1 (SVOC, Pesticides, Ignitability, pH, Oil & Grease, Paint Filter)
Sample Type:	FB (FIELD BLANK)

At the bottom of the form is a "Save Task" button.

Figure 52

There are three sorting options for the data displayed in the edit screen, indicated by the buttons at the top of the page. The functions of sort match the text on the button. For example, the 'Sort by Location' button orders the display of the task map by ascending location identifier. All data entry on the task map edit screen is not saved until the user clicks on the 'Click Here to Continue' icon.

When the user clicks on this icon, some validations are performed on the entered data. If any problems are found, the user is informed via a pop-up alert. Otherwise, another screen is displayed that will confirm that the records have been saved. This screen is shown in Figure 53.

In addition to the status message, two more buttons are added to the top of the screen. The second button on this screen, 'Return to the FCTM Main Screen', will return the user to the main screen of the Field Task Map application. The first, 'Ready to Load Data for this FCTM', will attempt to perform the loading function on the currently displayed FCTM. This is equivalent to selecting the 'Load Data' radio button and picking a field task map from the Field Task Map main screen (Figure 49). Clicking on this button will execute a series of checks on the data to determine if it is valid for the PEMS system.




Figure 53

These checks include, but are not limited to, a check for the uniqueness of the sample id, for the uniqueness of the form number and for the mandatory sample associations. An example of the output of this load is displayed in Figure 54.

If the loading operation is successful, the field task map is removed from the drop-down list on the Field Task Map main screen.




Figure 54

4.2.3 Field Chain of Custody (FCOC)

Clicking on the 'Click Here to Continue' icon on the loading status screen will display the main screen of the Field Chain of Custody application (See Figure 55), the next step in the sample planning process. The FCOC application is displayed because the majority of the samples are placed on a Field Chain of Custody form instead of the Flow Paced Composite form.

The FCOC main screen is divided into editing and printing operations. FCOC documents are edited by selecting a FCOC document from the FCOC Numbers drop-down list and the 'Edit FCOC' radio button (the default for this screen) and clicking on the 'Click Here to Continue' icon. By default, the FCOC Numbers drop-down list defaults to 'Create A New Field Chain Of Custody'. Since the loading of the task map created a FCOC, that FCOC will be used to illustrate the editing capabilities of the FCOC application.




Figure 55

Figure 56 shows the initial edit screen of the Field Chain of Custody application. This screen is used to edit the information that will be common to all of the samples listed on the FCOC document. The Delete FCOC option will not be displayed if any of the samples on the FCOC are on a finalized LCOC or have results loaded. 'The Logbook No.' field is used to indicate the logbook used to record the sampling operation, if any. The Team Leader drop-down list corresponds to a list of personnel in the 'Badged Person' support table. If the team leader is not listed in the drop-down list, that person can be added by using the 'Badged Person' application (refer to the Global Administration section of this document for more information). The 'No Team Leader' selection is the default value for this field.

The screenshot shows a Netscape browser window titled "PEMS - Netscape" with the address bar displaying "https://www-pems.bechteljacobs.org/qpems/qpems/Pems.load_menus". The left sidebar contains a navigation menu with links: Instruction, Manual Project, Lab Data, Data, Management, Global Admin, Project Admin, Reports, What's New, Switch Project, and Help. The main content area is titled "Edit FCOC" and includes a note: "Required fields are marked with +". Below this is a form with the following fields:

Field Chain of Custody (FCOC) Number:	froc_doc1 +	Rename FCOC
Logbook No.:	12413	
Planned Date(mm/dd/yyyy):	01/28/2003 +	
Team Leader:	703102(Eckerman, Mike)	

At the bottom of the form is a button labeled "CLICK HERE TO CONTINUE".

Figure 56

There are 2 additional functions on this screen that are represented by links. The 'Rename FCOC' link allows the user to change the document number for the FCOC. Clicking on this link will display a popup prompt for the new FCOC document number, as shown in Figure 57. The new FCOC number cannot be a number that is used by any other FCOC for the current project. Clicking on the 'OK' button

will rename the FCOC to the new document number. Clicking on the 'Cancel' button will return the user to the initial edit screen.

The 'Delete FCOC' link will display a prompt asking the user to confirm the operation, as shown in Figure 58. Clicking on the 'OK' button will delete the FCOC and all of the samples associated with it from the PEMS system. Clicking on the 'Cancel' button will return the user to the initial edit screen.




Figure 57




Figure 58

Clicking on the 'Click Here to Continue' icon will display the Choose Sample screen for the Edit FCOC application (See Figure 59). The header information from the previous screen is displayed at the top of the page. The user interface consists of a drop-down list, a 'Back' button and a 'Click Here to Continue' icon. The 'Back' button will return the display to the previous screen. The 'Click Here to Continue' icon will execute the function currently displayed in the drop-down list. The drop-down list contains 2 command options, which are followed by a list of the samples currently associated with the FCOC. The

command options are ‘Create a New Sample ID’ (which is the default) and ‘Add Samples Using Task Map’.




Figure 59

When the ‘Create a New Sample ID’ option is selected and the icon is clicked, the Edit Sample screen is displayed, as shown in Figure 60. Depending on the resolution of the user’s PC, there may be a scroll bar on the right side of the screen. The same rules apply for sample id values as in previous applications. All of the fields displayed on this screen are fields that have been seen previously on the Field Task Reference and Field Task Map applications, with five exceptions. The first exception is the Clean Lot No. field. This field is used to hold the lot number used by the manufacturer of sample bottles to ensure that they are free from contaminants.

The next three of these exceptions are displayed on the screen as the Top Depth, Bottom Depth and Depth Unit fields. The Edit Sample screen for the FCOC application is the only place in PEMS where the sample depth information can be recorded. The Top and Bottom Depth fields are numeric and will

accept decimal values. The Depth Unit field is limited to either feet or inches. These fields are required for samples in a soil matrix taken at a borehole or geoprobe location

The final exception is the Comments field, which is used to enter additional information about the sample. Examples of what values this field might contain include, but are not limited to, special collection or handling instructions to the field samplers or reasons for the collection of the sample.

PEMS - Netscape

File Edit View Go Bookmarks Tools Window Help

https://www.pems.bechteljacobs.org/qpems/qpems/Pems.load_menus Go Search

Home My Netscape Search Bookmarks Free AOL & U... RealPlayer Ho... Welcome to U... Liquid Music N...

Instruction Manual Project

Lab Data Data Management Global Admin Project Admin Reports What's New Switch Project Help

Field Chain of Custody (FCOC) - Edit Sample

Required fields are marked with *

Save Data Save Data and Add Another Sample Back

Field Chain of Custody (FCOC) Number	fcoc_doc1
Sample ID	<input type="text"/>
Matrix	AA-Ambient Air
Task/Monitoring Location	format: task code/loc id (loc desc) - site
Sample Type	FB (FIELD BLANK)
Sampling Procedure/Revision	---/Rev. 0 (No sampling procedure)
Container Volume	<input type="text"/> Search
Container Unit	<input type="text"/>
Container Type	<input type="text"/>
Clean Lot No.	<input type="text"/>
Assoc. Sample ID	<input type="text"/> Search
Top Depth	<input type="text"/>
Bottom Depth	<input type="text"/>
Depth Unit	No depth unit
Collection Device	20-2020 Photoionization air monitor-PID
Collection Method	<input type="text"/>
Averaging Period	<input type="text"/>

Done

Figure 60

The Associate Sample field is the same as the displayed in the Field Task Map. However, a search button has been added. Clicking on this button will display a pop-up window that lists all of the regular samples that are planned for that location on the FCOC. The Associate Sample is only required when the sample type is Field Replicate ('FR')

The buttons on the Edit Sample screen work in the same way as those on the Field Task Reference edit screen.

The second command option, Add Samples Using Task Map, displays an intermediate screen that contains a drop-down list and a 'Click Here to Continue' icon as shown in Figure 61. The drop-down list contains all of the Field Task Maps for the current project, whether or not they have been successfully loaded.




Figure 61

After selecting the appropriate Task Map and clicking on the icon, the ‘Choose Task Reference Samples’ screen is displayed (See Figure 62). This screen contains all of the general fields of the selected Task Map and a check box at the left side of the screen. The “general” fields are those fields that do not uniquely identify the sample. Fields that are not displayed on this screen are the FCOC form number (specified by the FCOC that is being edited), sample id, associate sample id and sampling procedure. These fields may be entered on the Edit Sample screen. By default, all of the samples on the task map are pre-selected. Clicking on the checkbox will exclude the sample from the insertion process.

The screenshot shows a Netscape browser window with the URL https://www-pems.bechteljacobs.org/qpems/qpems/Pems.load_menus. The page title is "Choosing Task Reference Samples".

At the top, it displays: Task: 1407J Matrix: WD Date: 01/28/2003

On the left, there is a vertical menu with the following items: Instruction, Manual Project, Lab Data, Data, Management, Global Admin, Project Admin, Reports, What's New, Switch Project, and Help.

The main content area contains a table with the following data:

<input checked="" type="checkbox"/> Use this record	Monitoring Location	8002	Sample Type:	REG
	Site	K-25	Preservative:	---
	Collection Device	BAILER	Process:	
	Analysis Group	1227CHEM1	Averaging Period:	
	Collection Method			
	Container Volume/Unit/Type	10/OZ/POLY		

Below the table, there is a red button that says "CLICK HERE TO CONTINUE".

Figure 62

When all of the samples have been selected from the Task Map, click on the ‘Click Here to Continue’ icon. This will display ‘Adding Samples By Task Reference’ screen where the user can enter the sample id(s) for the record(s) that had been selected (See Figure 63). The identifying information for the Field Task Map is no longer displayed and the checkbox has been replaced with an input field for the sample id. After entering the sample id, click on the continue icon to complete the insertion of the new sample(s) into the FCOC and return to the main screen of the Field Chain of Custody application. Note

that the entered sample id will be converted to uppercase text by the insertion process. Also, even though the sample procedure in the Field Task Map is not displayed on the ‘Choosing Task Reference Samples’ screen, the procedure as defined in the Field Task Map is inserted with the sample information into the FCOC. The user may need to edit the inserted samples on the Edit Samples screen to add additional information such as sample depth, comments or associated sample ids.

Adding Samples By Task Reference

Sample ID:	Monitoring Location:	8002	Sample Type:	REG
TSKREFADD1	Site:	K-25	Preservative:	---
	Collection Device:	BAILER		
	Analysis Group:	1227CHEM1 (SVOC, Pesticides, Ignitability, pH, Oil & Grease, Paint Filter)		
	Collection Method:		Averaging Period:	
	Container Volume/Unit/Type:	10/OZ/POLY		

CLICK HERE TO CONTINUE

Figure 63

4.2.3.1 Printing FCOCs and Labels

Multiple FCOCs may be printed by selecting the starting and ending form number, clicking on the ‘Print FCOC’ radio button and then clicking on the continue icon. This will bring up the same save/open prompt as the print operation for a single FCOC. The ‘save’ operation will proceed as it did for printing a single document. The ‘open’ operation will open all of the selected documents with each FCOC being displayed on a different page in the overall file.

The remaining print operations for the FCOC application deal with the printing of barcode labels. To successfully print barcode labels, the user must have the barcode fonts installed on their PC. These fonts

are available from the PEMS Home Page, if they are not installed. From the PEMS Home Page, select the 'Barcode Downloads' option from the left pane, and then click on the 'Download Code 39 Bar Code Fonts' link. This will transfer a self-extracting ZIP file that contains the fonts to the user's PC. When the ZIP file is executed, the font files are created in the directory where the ZIP file was transferred. Installation of the fonts will vary depending on the operating system. For Windows NT/2000, these font files should be placed in the %systemroot%\Winnt\Fonts directory, which is typically C:\Winnt\Fonts. If the fonts are properly installed, any text written in the Code 39 font in a Word or similar word processor will appear as a barcode.




Figure 64

Now that an FCOC has been defined, it can be printed for distribution to the field samplers. To accomplish this, select the range of FCOC documents to print using the drop-down lists, click on the 'Print FCOC' radio button on the FCOC main menu and click on the continue icon. Each drop-down list contains all of the FCOC documents that exist for the current project, sorted alphanumerically by form number. Depending upon the browser being used to access PEMS, a prompt will be displayed to the user asking if the file should be saved or opened. This prompt will be similar to Figure 64, which is the prompt displayed for a Netscape version 4.7 browser. If the 'save' option is selected, PEMS will display another prompt asking for the name of the file to be created. The destination of the file is dependent upon the browser being used, but the file itself will be in Rich-Text format. If the 'open' option is selected, PEMS will use whichever application has been defined in the user's browser for viewing/printing Rich-Text formatted files. If the number of samples on the FCOC is large, the output will flow onto other pages. An example of a FCOC document as displayed by the word processor is shown as Figure 65.

BECHTEL JACOBS
Bechtel Jacobs Company LLC

Project Environmental Measurement System (PEMS)
Field Chain of Custody

Page: 1 for fcoc_doc1
Date: 03/03/2003 3

Project:	INSTMAN	Field Chain of Custody No:	fcoc_doc1	Planned Date:	28-JAN-03
Task:	1407J	Field Matrix:	WD		
Logbook:		Team Leader:			

Sample ID	Monit Loc/ Smpl Type	Sampling Proc./Rev	Smpl Container Vol/Type/Qty	Clean Lot No/ Assoc Smpl ID	Sample Date Time/Status	Depth Top/Bot/Unit	Equipment ID/ Sampling Dev	Presv/Desc Process	Sample Com Parameter An
test-tb1	8002 TB	ESP-301-1 0	1 ea		:P	/	DIP	None /None /	8081 pesticides
testamp11	8002 REG	ER/WAG2- SOP3202 0	10 OZ POLY 1 ea		:P	/	BAI	None /None /	SVOC, Pesticide pH, Oil & Greas
TSKREFAD D1	8002 REG	ER/WAG2- SOP3202 0	10 OZ POLY 1 ea		:P	/	BAI	None /None /	SVOC, Pesticide pH, Oil & Greas

Relinquished by: _____ Received by: _____

Page 1 Sec 1 1/1 At 2.5" Ln 1 Col 1 REC TRK EXT OVR English (U.S)

Figure 65

Users may print large barcode labels with the 'Print Large Barcodes' radio button or small barcode labels with the 'Print Small Barcodes' radio button. (See section 5.3) The large barcode labels are 3 inches high by 4 inches wide in size. The small barcode labels are 1 inch high by 2 inches wide. Selecting the large barcode radio button will display a selection screen, as shown in Figure 66. This screen lists all of the samples currently on the FCOC on the left and a diagram showing the number and position of labels that will fit on an 8 ½ x 11 page of label stock. With the appropriate label stock and the barcode fonts installed, labels can be printed using a standard laser printer, avoiding the need for a barcode printer. To make the most efficient use of label stock, click on the starting position of the first label in the diagram. By selecting the 'Open' option on this prompt, the browser will open the default application for Rich-Text files to display the barcode output. An example of this output is shown as Figure 67. The barcodes may be printed immediately from the word processor or the file can be saved for later printing, as needed.

Note that there may be multiple copies of a barcode label for a sample or samples on the output. This is not an error. During the Project Definition phase, Parameter Analysis Groups were selected for the project. At that time, the number of labels that should be printed for the group code was entered. The barcode printing program uses the Parameter Analysis Group code and the number defined by the project to determine how many labels should be printed for the sample.

Small barcode labels are printed in the same manner as the large barcode label. They, too, will print multiple copies depending on the definition of the Parameter Analysis Group for the samples being printed. On selecting this option, a selection screen will be displayed, with the exception that more labels may be printed on a page. Depending on project operations, barcode labels may or may not need to be printed.







Figure 66

VV5U4A09.RTF - Microsoft Word

File Edit View Insert Format Tools Table Window Help

114%

Avery Wizard Times New Roman 10 B I U

	<p>Instruction Manual Project</p> <p>Sample ID: TSKREFADD1</p>  <p>FCOC No.: fcoc_doc1 Matrix Type: WD Smpl Type: REG Analysis: SWOC, Pesticides, Ignitability, pH, Oil & Grease, Container: POLY 10 OZ (1 ea) Holding Time: N/A Preservative: --- Filtered: No Monitoring Location: 8002 Sample Date: _____ Sample Time: _____ Collected By: _____</p>
<p>Instruction Manual Project</p> <p>Sample ID: test-th1</p>  <p>FCOC No.: fcoc_doc1 Matrix Type: WD Smpl Type: TB Analysis:</p>	<p>Instruction Manual Project</p> <p>Sample ID: testsmpl1</p>  <p>FCOC No.: fcoc_doc1 Matrix Type: WD Smpl Type: REG Analysis:</p>

Page 1 Sec 1 1/1 At 0.7" Ln 1 Col 1 REC TRK EXT OVR English (U.S)

Figure 67

4.2.4 Flow Paced Monthly Composite (FMPC)

When the Field Task Map was created earlier in this manual, the user had the option to specify either a 'FAS-01 / FCOC' or 'SW-08 / Composite' form type. Usually, a FCOC form type will be specified, but FMPC operations are another possible sample plan. By selecting the FMPC form type, the samples and other planning information will be associated with an FMPC form.

To complete the FMPC form information, the user should select the 'Flow Paced Monthly Composite' application from the Data Management menu. This will display the main screen of the application, as shown below.




Figure 68

The drop-down list of the FPMC application shows all of the current FPMC forms available for editing. By default, the 'Edit FPMC' radio button is selected. Clicking on the 'Click Here to Continue' icon will display the Edit Screen for the FPMC form header information, as shown in Figure 69. After entering the header information and clicking on the continue icon, the detail edit screen will be displayed. The content of the detail edit screen may vary, depending upon whether or not the sample was specified as Preserved or Non-Preserved on the header screen. The screen for preserved samples is shown in Figure 70 and the screen for non-preserved samples is displayed as Figure 71.

PEMS - Netscape

File Edit View Go Communicator Help

Back Forward Reload Home Search Netscape Print Security Shop Stop

Bookmarks Location: https://www.pems.bechteljacobs.org/qpems/qpems/Pems.load_menus

Instant Message Internet EIM Sites Lookup New&Cool LMES Internal H Oracle Document CNET.com - Welc WebWHOS: Search Seagull's Therm

Flow-Paced Monthly Composite Header Info
Required fields are marked with *

FPMC Number:	fpmcform1 *
Task Code/Monitoring Location:	1407/J/26-005
Composite Date(mm/dd/yyyy):	
Composite Time(hhmm):	
Sampling Procedure/Rev.:	ER/WAG2-SOP3202 / Rev. 0 *
Logbook Number:	<input type="text"/> Search
Compositor Badge No:	No Team Leader
Weather:	<input type="text"/>
Comments:	<input type="text"/>
Preserved/Non-Preserved:	<input checked="" type="radio"/> Preserved <input type="radio"/> Non - Preserved

CLICK HERE TO CONTINUE

Instruction Manual Project

Lab Data

Data Management

Global Admin

Project Admin

Reports

What's New

Switch Project

Help

Document: Done

Figure 69

The two radio buttons in Figure 68 note ‘Preserved’ and ‘Non-preserved’ samples. This information can be stored in PEMS by clicking on the appropriate choice for the composite sample. Additional samples may be entered or deleted using this function.

PEMS - Netscape

File Edit View Go Communicator Help

Back Forward Reload Home Search Netscape Print Security Shop Stop

Bookmarks Location: https://www-pems.bechteljacobs.org/qpems/qpems/Pems.load_menus What's Related

Instant Message Internet EIM Sites Lookup New&Cool LMES Internal H Oracle Document CNET.com - Welc WebWHOS: Search Seagull's Therm

Preserved Sample Compositing

(A) Total Monthly Composite Volume of Preserved Samples (ml):

Composite Sample ID: Required

Delete Samples	Date Collected	Sample ID	(B) Flow Volume	Flow Units	(C) Percent of Total (B)/(D)	Aliquot Amount (ml) (A)*(C)
<div> <input type="button" value="Delete Composite Samples"/> <input type="button" value="Add Composite Samples"/> </div> <p>CLICK HERE TO CONTINUE</p>						

Instruction Manual Project

Lab Data

Data Management

Global Admin

Project Admin

Reports

What's New

Switch Project

Help

Document: Done

Figure 70




Figure 71

If composite samples will be split to go to different laboratories, that information may be entered using the “Adding Preserved Split Samples’ function as shown in Figures 72 and 73.




Figure 72

PEMS - Netscape

File Edit View Go Communicator Help

Back Forward Reload Home Search Netscape Print Security Shop Stop

Bookmarks Location: https://www.pems.bechteljacobs.org/gpems/gpems/Pems.load_menus What's Related

Add Split Sample
Required fields are marked with *

FPMC Number:	fpmcform1
Sample ID:	<input type="text"/> *
Monitoring Location:	26-005(EAST TENNESSEE TECHNOLOGY PARK) ▾
Sampling Dev.:	20-2020 Photoionization air monitor-PID ▾
Procedure No.:	ER/WAG2-SOP3202
Procedure Rev.:	0
Sample Type:	FB (FIELD BLANK) ▾
Collection Method:	<input type="text"/> ▾
Averaging Period:	<input type="text"/> ▾
Preserved:	--- (NO PRESERVATIVE REQUIRED) ▾
Analysis Group:	1227CHEM1 (SVOC, Pesticides, Ignitability, pH, Oil & Grease, Paint Filter) ▾
Container Volume:	<input type="text"/> Search
Container Unit:	<input type="text"/>
Container Type:	<input type="text"/>

CLICK HERE TO CONTINUE

Instruction Manual Project

Lab Data

Data Management

Global Admin

Project Admin

Reports

What's New

Switch Project

[Help](#)

Document: Done

Figure 73

After the information has been entered, use the 'Click Here to Continue' button to save the information. The screen depicted in Figure 74 will confirm that the composite information has been successfully saved to the database.




Figure 74

Figure 75 displays the Preserved BLANK FPMC Report sample query screen. By clicking on the 'Pick Samples' box, samples may be added to a FCOC.




Figure 75


Figure 76 shows FPMC BLANK Report output that may be printed for the samplers to use in the field.

fpmc_blank.RTF - Microsoft Word

File Edit View Insert Format Tools Table Window Help

Normal Times New Roman 10 B I U

2 3 4 5 6 7 8 9 10



**BECHTEL
JACOBS**
Bechtel Jacobs Company LLC

Project Environmental Measurement System (PEMS)
Flow-Paced Monthly Composite Form
Preserved Sample Compositing for January , 2003

Page: 1 of 1
Date: 03/21/2003 11:26 AM

Project: INSTMAN	Task: 1407J	Monitoring Loc: 26-005	
FPMC No.: fpmcform1			

Compositor:	FOP Number / Revision: ER/WAG2-SOP3202/0	Date:	
Field Logbook:	Weather: Cloudy	Time:	

(A) Total Monthly Composite Volume of Preserved Samples (mL):

Weekly Retrieval Number	FCOC Number	Planned Date	Date Collected	Sample ID	(B) Total Flow Volume Since Last Retrieval (gal)	(C) Percent of Total Flow Volume (B/D)*100	Aliquot Amount (mL) A*C/100
4	fcoc_doc1	01/28/2003		TSKREFAD D1			
4	fcoc_doc1	01/28/2003		test-th1			
4	fcoc_doc1	01/28/2003		testsmpl1			

(D) Total Flow Volume for Month (gal): 0

Comments:

Composite Sample ID:

Surface Water Split Sample Compositing

Sample ID	Sample Time	Associated Sample ID	Sample Type	Sample Status	Container Volume	Container Type	Parameter Analysis Group	pH	Preservative Code

Page 1 Sec 1 1/1 At 1.5" Ln 1 Col 1 REC TRK EXT OVR English (U.S.)

Figure 76

5.0 Field Data

5.1 Enter Collection Date and Time

Depending on the operations of the project, samples may be planned days, weeks or even months in advance of the actual sample collection operations. When the samples are finally collected, the collection date and time must be recorded in PEMS. This is accomplished by using the 'Enter Collection Date And Time' application from the Data Management menu. The main screen of this application is show below, as Figure 78.

The screenshot shows a Netscape browser window titled 'PEMS - Netscape'. The address bar displays 'https://www.pems.bechteljacobs.org/gpems/gpems/Pems.load_menus'. A left-hand navigation menu contains links: 'Instruction Manual Project', 'Lab Data', 'Data Management' (highlighted in yellow), 'Global Admin', 'Project Admin', 'Reports', 'What's New', 'Switch Project', and 'Help'. The main content area is titled 'Enter Collection Date and Time'. It features a form with a green header 'Field Chain of Custody (FCOC) Numbers' containing a dropdown menu with 'fcoc.doc1' selected. Below this is a green box with the text 'Enter Sample Id (Optional)' and 'Leave it blank if all Samples are desired', followed by a white input field. A note states: 'For FCOCs that contain many samples, the above Sample field can be used to reduce the Sample selection list. This field is not case sensitive and partial Sample Ids may be entered.' Another note below says: 'Note: FCOC form numbers will not be displayed if all samples on the FCOC are marked as lost, retired, uncollectible or if all samples are on a finalized LCOC.' At the bottom of the form is a red button with white text that says 'CLICK HERE TO CONTINUE'.

Figure 78

Users are able to enter date/time collected for one (by entering the sample id - see Figure 79) or multiple samples (by leaving the sample id field blank - see Figure 80). After entering the date and time and selecting the sample, clicking on the 'continue' icon, will save the information for the sample.




Figure 79

To enter the data/time for multiple samples associated with a FCOC, select the samples to edit (see Figure 80). This screen lists all of the samples for the selected FCOC with checkboxes to include samples in the collection update. After the samples have been specified and the collection data entered, clicking on the 'continue' icon will save the information for the sample.

Using the Collection Date And Time application will also change the status of the sample. When the collection information is saved, the sample status is changed 'C'

(Collected). The result of changing the sample status is that that sample may now be placed on a Laboratory Chain of Custody (LCOC).




Figure 80

5.2 Laboratory Chain of Custody (LCOC)

LCOCs are used to document the transfer of collected samples to the laboratory for analysis. To execute this application, click on the Laboratory Chain of Custody link on the Data Management menu. The main screen of the LCOC application is displayed in Figure 81. This screen contains a drop-down list that defaults to the 'Create A New Lab Chain of Custody (LCOC)'. This drop-down list also shows all of the LCOCs by form number for the current project, the SOW that is assigned to the LCOC and the finalized status of the LCOC. The finalized status is either 'Finalized' or null (indicating that the LCOC is not finalized).

Some projects prepare their LCOC documents in advance of the actual collection of samples to more rapidly facilitate shipping. In order to do this, PEMS must store the LCOC and allow editing of the information to reflect actual operations. Finalizing an LCOC locks the document to prevent editing and indicates to the system that the samples have been shipped for analysis. Analytical results may only be loaded for samples that are on finalized LCOCs.




Figure 81

To create/edit an LCOC, ensure that the Edit LCOC radio button is selected and click on the 'Click Here to Continue' icon. If the LCOC is being created, the LCOC form screen will be displayed as shown in Figure 82. If the LCOC is being edited, a link will appear next to the Form Number field. This link, labeled 'Rename this LCOC', will display a prompt asking the user for the new form number. Entering the new form number and clicking on the 'OK' button will rename the LCOC to the new form number.

Another difference in the edit screen from the create screen is the addition of the 'Delete LCOC' button at the bottom of the screen. Clicking on that button will delete the LCOC and mark the samples currently on that LCOC as available for placement on a new LCOC.

The user should enter the LCOC information into the fields present. Note that the finalization check box beneath the form number field. By clicking on this check box, the LCOC will be finalized when the application is completed.

Edit LCOC

Required fields are marked with *

LCOC Form Number:	<input type="text"/> *
	If finalized <input type="checkbox"/>
Lab Ship Date(mm/dd/yyyy):	<input type="text"/>
Lab:	<input type="text"/> * <input type="button" value="Search"/>
SOW:	<input type="text"/> *
SOW Rev.:	<input type="text"/> *
Ship Container No.:	<input type="text"/> *
Cooler Temp (C):	<input type="text"/>
Airbill Number:	<input type="text"/>

CLICK HERE TO CONTINUE

Figure 82

The Shipment Group identifier defaults to the form number and is displayed on the printed LCOC. The Search button displayed next to the Lab field will access all of the SOWs defined for the current project and display them in a pop-up window, as shown in Figure 83. It is through this window that the samples are assigned to the laboratories that will perform the analysis.

Laboratory/ SOW /SOW Rev.

DCSL - BKETO1320 - 0

OK Cancel

Figure 83

After entering all of the form information and clicking on the continue icon, the LCOC Choose Samples screen is displayed as shown in Figure 84.




Figure 84

After selecting the samples to be placed on the LCOC by clicking on the checkboxes next to the sample id, clicking on the continue icon will add the samples to the LCOC and return the user to the main screen of the application.

5.2.1 Printing LCOCs

Using the 'Print LCOC' function, the LCOC will be in a printable Word format as shown in Figure 85. LCOCs must be printed one at a time.

lcoc.RTF - Microsoft Word

File Edit View Insert Format Tools Table Window Help

Normal Times New Roman 16 B I U

1 2 3 4 5 6 7 8 9 10

Project Environment Measurements System (PEMS)

Laboratory Lab Chain of Custody (LCOC)
 Project: INSTMAN SOW: BKET01320
 Lab Chain of Custody (LCOC) Number: INST_1

Date: 03/21/2003 03:38 PM

Sample ID	FCOC	Task	Matrix	Smpl Type	Vol	Unit	Smpl Container Type	Qty	Sample Date	Smpl Time	Assoc Smpl ID	Param Anly Group	Preservative	Remarks	
TSKREFADD1	fcoc_docl	1407J	WD	REG	10	OZ	POLY	1	21-MAR-2003	1249		SVOC, Pesticides, Ignitability, pH, Oil & Grease, Paint Filter	None		
Laboratory:			Date Submitted To Lab:			Ship Container No:			SSG ID:		Cooler Temperature:		Airbill No:		SOW
DCSL - DCSL			21-MAR-2003			COOLER1			1		DEG C				BKET01320
Relinquished by:			Date(mm/dd/yyyy):			Time:			Received by:			Date(mm/dd/yyyy):		Time:	
Relinquished by:			Date(mm/dd/yyyy):			Time:			Received by:			Date(mm/dd/yyyy):		Time:	
Relinquished by:			Date(mm/dd/yyyy):			Time:			Received by:			Date(mm/dd/yyyy):		Time:	
Relinquished by:			Date(mm/dd/yyyy):			Time:			Received by:			Date(mm/dd/yyyy):		Time:	
Comments:															

—

Page 1 Sec 1 1/1 At 0.5" Ln 1 Col 1 REC TRK EXT OVR English (U.S)

Figure 85

Barcode sample bottle labels can also be printed from the LCOC menu.

5.2.2 Electronic Export LCOC

From the Laboratory Chain of Custody main screen, the 'Electronic Export LCOC' button will download a .csv file to the user's PC as shown in Figure 86.

Microsoft Excel - lcoc.csv

File Edit View Insert Format Tools Data Window Help

A1 = EP_ID

	A	B	C	D	E	F	G	H	I	J	K	L	M	N	O	P
1	EP_ID	LCOC_No	Smpl_ID	FCOC_No	Task	Matrix	Smpl Type	Smpl_Con	Smpl_Con	Smpl_Con	Smpl_Con	Smpl Date	Smpl Time	Assoc_Sm	Param_An	Pre
2	INSTMAN	INST_1	TSKREFA	fcoc_doc1	1407J	WD	REG	10	OZ	POLY	1	3/21/03	1249		SVOC, Pe...	
3																
4																
5																
6																
7																
8																
9																
10																
11																
12																
13																
14																
15																
16																
17																
18																
19																
20																
21																
22																
23																
24																
25																
26																
27																
28																
29																
30																
31																
32																
33																
34																
35																
36																
37																

PDF 4.0 Prompt

lcoc

Ready

NUM

Figure 86

5.3 Bar Code Operations

Note: Please contact PEMS support for help with barcode readers, programs, etc. The PEMS Bar Code User is posted on the PEMS website at the following location:
https://www-pems.ettp.energy.gov/pems/home/PEMS_Bar_Code_Users_Guide.pdf

PEMS Support Email: pems@ettp.doe.gov

Download PEMS.HEX from the PEMS home page.
<https://www-pems.ettp.energy.gov/pems/home/home.html>

Users need an Internet-ready PC with a 9-pin male COM1 port to load C Programs into the bar code readers. The user will need to download the COMMUNICATIONS.EXE file from the PEMS homepage, then unzip the file and follow the readme.txt instructions.

The project may use barcode readers for collection of field data and sample data. During the course of the operation the Bar Code reader, two files will be generated:

FIELD.DAT, A file containing field measurement transactions.

SAMPLE.DAT, A file containing sample collection transactions.

Each transaction appends a row to one of the files, depending on the type of transaction. The files are in a column-oriented format (fixed position). This file is uploaded to your PC and then uploaded in the PEMS database. The file may be modified while it is on the user's PC before uploading it to the PEMS database. The file formats are described below.

FIELD.DAT

Beginning Column Position	Length	Field Name	Description
1	7	*Reader No	The Bar Code Reader No (i.e. PDT0100)
9	10	*Project ID	The Project ID
20	7	*Badge No	The Team Leader Badge No
28	15	*Logbook No	The Logbook No
44	2	*Matrix	The Matrix Type (see PEMS for valid values)
47	9	Equipment ID	The Equipment ID
57	15	*Monitoring Location	The Monitoring Location
73	15	*Form ID	The Field Form ID
89	5	*Task Code	The Task Code
95	10	*Planning Date	The Planning Date in (MM/DD/YYYY)
106	9	*Parameter	The Field Measurement Parameter Code (see PEMS for valid values)
116	16	*Field Units	The Field Measurement Unit
133	10	*Field Result	The Field Measurement Result in numeric decimal format
144	20	*Weather	The Weather
165	20	*Comments	The Comments
186	10	Field Date	The Field Measurement Collection Date in (MM/DD/YYYY)
197	4	*Field Time	The Field Measurement Collection Time in (HHMM)

* Required fields

SAMPLE.DAT

Beginning Column Position	Length	Field Name	Description
1	7	*Reader No	The Bar Code Reader No (i.e. PDT0100)
9	10	*Project ID	The Project ID
20	7	*Badge No	The Team Leader Badge No
28	15	*Logbook No	The Logbook No
44	2	*Matrix	The Matrix Type (see PEMS for valid values)
47	9	Equipment ID	The Equipment ID
57	15	*Monitoring Location	The Monitoring Location
73	15	*Sample ID	The Sample ID (for each container)
89	15	Lot No	The Certified Lot No
105	10	Top Depth	The Top Depth in numeric decimal format
116	10	Bottom Depth	The Bottom Depth in numeric decimal format
127	16	Depth Units	The Depth Units
144	10	*Sample Date	The Sample Date in (MM/DD/YYYY)
155	4	*Sample Time	The Sample Time in (HHMM)

* Required fields

5.3.1 Print Bar Code Menu Sheets

The data may be uploaded into the database from a personal computer. Barcodes for field measurements are displayed below.




Figure 87

Figure 88 displays some of the barcodes available in the system for field measurements.




Figure 88

The user can print barcodes by parameter in the order desired. The user should enter the order in the fields as shown in Figure 89.




Figure 89

Barcodes may also be printed to enter field measurements into PEMS as shown in Figure 90.




Figure 90

5.3.2 Transfer Barcode Reader Data

The data entered into the barcode reader is downloaded to a PC, then uploaded into PEMS. The file name may be entered into the field or the browser may be used to pick the file name. The file format must match the radio button selected (Figure 91).




Figure 91

The system runs checks on the data so that the user can be notified immediately if there are errors in the data (Figure 92). The user must correct the file and re-transfer it to PEMS.




Figure 92

When the data is correct, it will be displayed as in Figure 93. The "Click Here to Continue" button will save the data to the database.

PEMS - Netscape

File Edit View Go Communicator Help

Back Forward Reload Home Search Netscape Print Security Shop Stop

Bookmarks Location: https://www.pems.bechteljacobs.org/qpems/qpems/Pems.load_menus What's Related

Instant Message Internet EIM Sites Lookup New&Cool LMES Internal H Oracle Document CNET.com - Welc WebWHOS: Search Seagull's Therm

Instruction Manual Project
Lab Data
Data
Management
Global Admin
Project Admin
Reports
What's New
Switch Project
Help

Loading Barcode Sample Data

File ID: B0016620
Original File name: CAWorking\B0039412.html

File Name	Database Name	Line No	Reader No	Project ID	Badge No	Logbook No	Matrix Type	Equipment ID	Monitoring Location	Sample ID	Cleaning Lot No	Top Depth	Bottom Depth	Depth Units
B0039412.HTML	B0016620	1	PDT0101	INSTMAN	703102	12413	WS		WB 9619-5	WB-01				
B0039412.HTML	B0016620	2	PDT0101	INSTMAN	703102	12413	WS		WB 9619-5	WB-02				
B0039412.HTML	B0016620	3	PDT0101	INSTMAN	703102	12413	WS		WB 9619-5	WB-03				
B0039412.HTML	B0016620	4	PDT0101	INSTMAN	703102	12413	WS		WB 9619-5	WB-04				
B0039412.HTML	B0016620	5	PDT0101	INSTMAN	703102	12413	WS		WB 9619-5	WB-05				
B0039412.HTML	B0016620	6	PDT0101	INSTMAN	703102	12413	WS		WB 9619-5	WB-06				
B0039412.HTML	B0016620	7	PDT0101	INSTMAN	703102	12413	WS		WB 9619-5	WB-07				
B0039412.HTML	B0016620	8	PDT0101	INSTMAN	703102	12413	WS		WB 9619-5	WB-08				
B0039412.HTML	B0016620	9	PDT0101	INSTMAN	703102	12413	WS		WB 9619-5	WB-10				
B0039412.HTML	B0016620	10	PDT0101	INSTMAN	703102	12413	WS		WB 9619-5	WB-09				
B0039412.HTML	B0016620	11	PDT0101	INSTMAN	703102	12413	WS		WB 9619-5	WB-11				
B0039412.HTML	B0016620	12	PDT0101	INSTMAN	703102	12413	WS		WB 9619-5	WB-12				
B0039412.HTML	B0016620	13	PDT0101	INSTMAN	703102	12413	WS		WB 9619-5	WB-13				
B0039412.HTML	B0016620	14	PDT0101	INSTMAN	703102	12413	SO		WB 9619-5	WB-14				
B0039412.HTML	B0016620	15	PDT0101	INSTMAN	703102	12413	SO		WB 9619-5	WB-15				
B0039412.HTML	B0016620	16	PDT0101	INSTMAN	703102	12413	SO		WB 9619-5	WB-16				

CLICK HERE TO CONTINUE

Figure 93

5.3.3 Edit Sample Barcode Reader Data

Barcode data may be edited by using the Edit function and picking the FCOC number to edit.




Figure 94

The FCOC data will be displayed. (Figure 95) Each row to be edited may be chosen separately.




Figure 95

The data may be edited or deleted as needed as shown in Figure 96.

PEMS - Netscape

File Edit View Go Communicator Help

Back Forward Reload Home Search Netscape Print Security Shop Stop

Location: https://www.pems.bechteljacobs.org/qpems/qpems/Pems.load_menus

Instant Message Internet EIM Sites Lookup New&Cool LMES Internal H Oracle Document CNET.com - Welc WebWHOS: Search Seagull's Therm

Instruction
Manual Project

Lab Data
Data
Management
Global Admin
Project Admin
Reports
What's New
Switch Project
Help

Edit Sample Barcode Reader Data

Project ID	Filename	File ID	Line
INSTMAN	B0039412.HTML	B0016620	1

Delete Record

Required fields are marked with *

Reader No:	PDT0101 *
Team Leader:	703102(Eckerman, Mike) *
Logbook No:	12413 *
Matrix Type:	WS
Equipment ID:	None
Monitoring Location:	26-005
Sample ID:	WB-01 [ML= FCOC=] *
Cleaning Lot:	
Top Depth:	
Bottom Depth:	
Depth Units:	None
Sample Date(mm/dd/yyyy):	01/16/2003 *
Sample Time(HHMM):	0930 *

**CLICK HERE
TO CONTINUE**

Figure 96

The system will display the corrections that should be made to the data before it is loaded into the database. (Figure 97)

PEMS - Netscape

File Edit View Go Communicator Help

Back Forward Reload Home Search Netscape Print Security Shop Stop

Bookmarks Location: https://www.pems.bechteljacobs.org/qpems/qpems/Pems.load_menus

Instant Message Internet EIM Sites Lookup New&Cool LMES Internal H Oracle Document CNET.com - Welc WebWHOS: Search Seagull's Therm

Instruction Manual Project

Lab Data

Data Management

Global Admin

Project Admin

Reports

What's New

Switch Project

Help

Fix the fields indicated below.

Project ID	Filename	File ID	Line
INSTMAN	B0039412.HTML	B0016620	1

Delete Record

Required fields are marked with *

Reader No:	PDT0101 *
Team Leader:	703102(Eckerman, Mike) *
Logbook No:	12413 *
Matrix Type:	WS
Equipment ID:	None
Monitoring Location:	26-005
Sample ID:	WB-01 [ML= FCOC=] *
↑ A sample must be planned before applying barcode results	
Cleaning Lot:	
Top Depth:	
Bottom Depth:	
Depth Units:	None
Sample Date(mm/dd/yyyy):	01/16/2003 *
Sample Time(HHMM):	0930 *

CLICK HERE TO CONTINUE

Document: Done

Figure 97

When the barcode data is complete and correct, the system will indicate that it has been saved to the PEMS database as shown in Figure 98.




Figure 98

5.4 Field Measurements

From the Data Management main screen, the user may click on Field Measurements. The three radio buttons allow the user to enter any field measurements taken during sampling activities, print them or delete as appropriate. (Figure 99)




Figure 99

Users have several options for field measurement data entry.

1. Create a new field measurement by directly entering field data.
2. Create a new template. This feature is good for on-going projects collecting a standard set of field measurements.
3. Managing existing templates

Templates are easy organization tools if the project has standard sets of field measurement. They can be set up in any order and can be renamed. It is especially

helpful in ensuring the units are correct per project requirements. Templates can be created for the current project only. (See Figure 100) Name the template and pick the units associated with that field measurement. They can be changed by hand during data entry. Using the 'Click Here to Continue' saves the data.

PEMS - Netscape

File Edit View Go Communicator Help

Back Forward Reload Home Search Netscape Print Security Shop Stop

Bookmarks Location: https://www-pems.bchteljacobs.org/qpems/qpems/Pems.load_menus What's Related

Instant Message Internet EIM Sites Lookup New&Cool LMES Internal H Oracle Document CNET.com - Welc WebWHOS: Search Seagull's Therm

Instruction
Manual Project

Lab Data
Data
Management
Global Admin
Project Admin
Reports
What's New
Switch Project
Help

Adding Template

**CLICK HERE
TO CONTINUE**

Template Name:	
Parameter	Unit
<input type="text"/>	<input type="text"/> Search
<input type="text"/>	<input type="text"/> Search
<input type="text"/>	<input type="text"/> Search
<input type="text"/>	<input type="text"/> Search
<input type="text"/>	<input type="text"/> Search
<input type="text"/>	<input type="text"/> Search
<input type="text"/>	<input type="text"/> Search
<input type="text"/>	<input type="text"/> Search
<input type="text"/>	<input type="text"/> Search
<input type="text"/>	<input type="text"/> Search
<input type="text"/>	<input type="text"/> Search

**CLICK HERE
TO CONTINUE**

Document: Done

Figure 100

Templates can be copied from one project to another project by using the copy function and making the new template name unique. The user must be logged into the project to which the template will be copied.

A template can be deleted if it doesn't suit the user's current needs or a mistake is made in the entry process.




Figure 101

The user may choose the appropriate template from a drop down list of templates for a project.




Figure 102

The 'Edit Generic Field Measurement' screen allows the user to enter generic field measurements by entering header information applicable to one or multiple individual field measurements. (Figure 103)

Edit Generic Field Measurement

Required fields are marked with ➔

CLICK HERE TO CONTINUE

Form ID:	<input type="text"/> ➔
Configuration Control:	No ➔
Monitoring Location/Task:	26-005 / 1407J ➔
Begin Date(mmm/dd/yyyy):	<input type="text"/> ➔
Field Logbook Number:	<input type="text"/> Search
Matrix Type:	AA - Ambient Air ➔
Team Leader:	No Team Leader
Collection Date(mmm/dd/yyyy):	<input type="text"/> ➔
Collection Time(HHMM):	<input type="text"/> ➔
Weather:	<input type="text"/>
Comments:	<input type="text"/>

CLICK HERE TO CONTINUE

Figure 103

Only the fields with arrows are mandatory. The form ID number must be unique within the project. The configuration control option allows the field measurement to be “locked” to prevent changes. This option should be used only after data entry and review. The begin date and collection date are frequently the same (i.e., the date measurement was taken in the field). The team leader field has a pick list of names that are populated in the Global Admin/Badged Person section of PEMS.

Select ‘Add New Field Parameter(s)’ to enter field measurement results.




Figure 104

To directly enter individual field measurements into the PEMS database, users can select the parameter name(s) from the drop down list and enter the results and appropriate units. (Figure 105)

PEMS - Netscape

File Edit View Go Communicator Help

Back Forward Reload Home Search Netscape Print Security Shop Stop

Bookmarks Location: https://www.pems.bechteljacobs.org/qpems/qpems/Pems.load_menus What's Related

Instant Message Internet EIM Sites Lookup New&Cool LMES Internal H Oracle Document CNET.com - Welc WebWHOS: Search Seagull's Therm

Instruction Manual Project

Lab Data

Data Management

Global Admin

Project Admin

Reports

What's New

Switch Project

Help

Add Generic Field Measurements

CLICK HERE TO CONTINUE

Parameter	Result	Unit	Time	Equipment ID
Conductivity		OHM	1557	
Depth		FT	1557	
			1557	
			1557	
			1557	
			1557	
			1557	
			1557	
			1557	
			1557	
			1557	

CLICK HERE TO CONTINUE

Figure 105

A single field measurement may be edited if needed as shown in Figure 106.




Figure 106

5.5 Associate Samples

Regular samples can be associated or linked with the appropriate quality control samples.
Select the desired QC sample id.




Figure 107

Select the 'Add Associate Sample ID(s) to this Sample' option as shown in Figure 108.




Figure 108

Click on the check boxes for regular samples to associate with the appropriate QC sample.




Figure 109

PEMS will display a table showing the associations the user created. Associations between regular and QC samples can be removed by using the delete function shown in Figure 110.




Figure 110

5.6 Associate Samples to Field Measurements

Samples can be associated to field measurements using this function. The association criteria are Location and Matrix. Select the desired sample ID and the appropriate criteria.

The screenshot shows a Netscape browser window titled "PEMS - Netscape". The address bar displays the URL "https://www.pems.bechteljacobs.org/qpems/qpems/Pems.load_menus". The browser's toolbar includes buttons for Back, Forward, Reload, Home, Search, Netscape, Print, Security, Shop, and Stop. Below the toolbar, there are links for Bookmarks, Location, and What's Related. The main content area features a blue sidebar on the left with a menu containing: Instruction Manual Project, Lab Data, Data Management (highlighted in yellow), Global Admin, Project Admin, Reports, What's New, Switch Project, and Help. The main content area has a title "Associate Samples To Field Measurements Query Screen" and a subtitle "Choose a Sample ID". Below the subtitle is a table with the following data:

Sample ID	Type	Task	Location	Anly Group	Date
TSKREFADD1	REG	1407J	8002	1227CHEM1	21-MAR-2003
fpmc1	REG	1407J	26-005	1227CHEM1	
split1	FB	1407J	26-005	1227CHEM1	
test-tb1	TB	1407J	8002	8081PEST	21-MAR-2003
testsmpl1	REG	1407J	8002	1227CHEM1	21-MAR-2003

Figure 111

Next, select the 'Add Associate Field Measurement to the Sample' option (Figure 112).




Figure 112

Check the box next to the applicable field measurement. Save the information by selecting 'Click here to continue'.




Figure 113

Associations between samples and field measurements can be removed using the delete option. Click in the box of the row to be deleted, then “Click Here to continue” to save the information.




Figure 114

Comments may be added through the ‘Revise comments to field Meas/Asoc Sample’ option shown in Figure 115.




Figure 115

6.0 Laboratory Data

All laboratories must use the AMSED electronic data deliverable (EDD) format in order to load data into PEMS. This is accomplished utilizing the application under the Lab Data menu in the left frame of the PEMS system. Users with the proper roles (e.g., PDC) can retrieve copies of EDD files using the “Get a Hard Copy of EDD Data From the Database” button and can view the file history and SOW information. These users can get a missing analytes report to ensure all results are loaded into PEMS.

6.1 Transfer and Validate EDD Data

The appropriate laboratory will log on to PEMS and load data files from this menu using the ‘Transfer and Validate EDD Data’ application shown in Figure 117.




Figure 117

The laboratory can enter the file name in the field or may use the ‘Browse’ button. When the file name and path are entered, click on the “Load the Data” button as shown in

Figure 118. The laboratories will use the PEMS file naming convention for all files as stated in the AMSED specification. Radiological file names must begin with an “r”, and end with the first seven characters of the batch id or SDG number. The five radiological file extensions are .res (results), .lcs (laboratory control sample, matrix spike, and matrix spike duplicates), .dup (duplicate), .tir (tentatively identified radionuclides), and .mb (blank file format). Environmental or ‘non-rad’ files must begin with an “n”, and end with the first seven characters of the batch or SDG number. The file extensions are .res (results), .lcs (laboratory control sample), .ms (matrix spike, matrix spike duplicate and duplicate), and .tic (tentatively identified compounds).




Figure 118

PEMS performs validation checks against the data supplied in the EDD file to ensure that it matches system-specific codes and uses valid values for numeric and date fields. There

are two error categories, Warning and Fatal. A warning message means that contents of a field are not exactly what is required in PEMS but will allow the data to load. A fatal message means that the data will not load until those errors are corrected. The EDD Loader program will let the laboratory user know the error and the row the error occurred as seen in Figure 119. The laboratory may use this feature as many times as it takes for the data to load.

Consistency Check Radiological Results Data

Database name: L0016621

Main Menu

File: C:\Working\r0023243.res

You have fatal errors that must be corrected

Record #	Field	Error	Error/Warning
1	TSKREFADD	No LCOC for Client Sample ID Search	Fatal
2	TSKREFADD	No LCOC for Client Sample ID Search	Fatal
2	Total Strontium	Nonstandard Analyte Name	Warning
3	TSKREFADD	No LCOC for Client Sample ID Search	Fatal
4	TSKREFADD	No LCOC for Client Sample ID Search	Fatal

Figure 119

If no fatal errors are encountered during the validation of the data file, a success screen is displayed as shown in Figure 120. The user should click on the 'Save Data To Database' button to save the EDD data to the PEMS database. If this button is not clicked at this point, then the data will not be saved.




Figure 120

When the data file is loaded into PEMS, a database name will be assigned to that file. The user may click on the “OK” button to leave this screen, or may view the data that was just loaded as shown in Figure 121.




Figure 121

6.2 Get a Hardcopy of EDD Data from the Database

A user can retrieve a copy of a project EDD file from this screen. A date range is provided to assist in locating the desired data file. The user can retrieve data from the current or all projects to which the user has access. See Figure 122. The default date range is two weeks and the value of the date range can be changed.

The screenshot shows a Netscape browser window with the title 'PEMS - Netscape'. The address bar displays the URL 'https://www.pems.bechteljacobs.org/qpems/qpems/Pems.load_menus'. The page content includes a blue sidebar on the left with the following links: 'Instruction Manual Project', 'Lab Data', 'Data Management', 'Global Admin', 'Project Admin', 'Reports', 'What's New', 'Switch Project', and 'Help'. The main content area has the heading 'Get A Hardcopy of EDD Data From The Database'. Below the heading is a form with three rows: 'Start Date(mm/dd/yyyy):' with a text box containing '03/07/2003', 'End Date(mm/dd/yyyy):' with a text box containing '03/21/2003', and 'Project:' with two radio button options, 'Current Project' (selected) and 'All Projects'. Below the form is a red button labeled 'CLICK HERE TO CONTINUE'.

Start Date(mm/dd/yyyy):	03/07/2003
End Date(mm/dd/yyyy):	03/21/2003
Project:	<input checked="" type="radio"/> Current Project <input type="radio"/> All Projects

CLICK HERE TO CONTINUE

Figure 122

Figure 123 shows a list of files displayed by using the query shown in Figure 122. The files can be downloaded to the user's PC by clicking on the 'Database Name' field. The application will allow the user to save the file on the PC at the user's discretion. To retrieve a file, click on the database name of the file. The PC window will appear asking for the directory in which to save the file. The filename can be changed at the user's discretion at this time. Clicking the "Save" button will save the file to the directory of choice.

PEMS - Netscape

File Edit View Go Communicator Help

Back Forward Reload Home Search Netscape Print Security Shop Stop

Bookmarks Location: https://www.pems.bechteljacobs.org/qpems/qpems/Pems.load_menus What's Related

Instant Message Internet EIM Sites Lookup New&Cool LMES Internal H Oracle Document CNET.com - Welc WebWHOS: Search Seagull's Therm

Instruction Manual Project

Lab Data

Data Management Global Admin Project Admin Reports What's New Switch Project Help

Get A Hardcopy of EDD Data From The Database Project: INSTMAN

Data Files				
Date Loaded	Database Name	Original File Name - Project - SOW Number(s).	File Type	Userid
Friday - 03/21/2003 05:47 pm	L0016624 Inserted	C:\Working\r0023243.res - INSTMAN - BKET01320	Rad Results	QA
Friday - 03/21/2003 05:45 pm	L0016623 Not Inserted	C:\Working\r0023243.res - INSTMAN	Rad Results	QA
Friday - 03/21/2003 05:43 pm	L0016622 Not Inserted	C:\Working\r0023243.res - INSTMAN	Rad Results	QA
Friday - 03/21/2003 05:31 pm	L0016621 Not Inserted	C:\Working\r0023243.res - INSTMAN	Rad Results	QA
Friday - 03/21/2003 04:52 pm	B0016620 Not Inserted	C:\Working\B0039412.html - INSTMAN	Sample	QA
Friday - 03/21/2003 04:47 pm	B0016619 Not Inserted	C:\Working\B0039412.html - INSTMAN	Sample	QA
Friday - 03/21/2003 04:46 pm	B0016618 Not Inserted	C:\Working\B0039412.html - INSTMAN	Sample	QA
Friday - 03/21/2003 04:44 pm	B0016617 Not Inserted	C:\Working\B0039412.html - INSTMAN	Sample	QA
Friday - 03/21/2003 04:42 pm	B0016616 Not Inserted	C:\Working\B0039412.html - INSTMAN	Field	QA

Document: Done

Figure 123

6.3 EDD Data File History

The user can retrieve a file by file name or by DRG number (also called "database name"), and see the history of that particular file if multiple versions have been submitted to the database to add or correct file data. The File History prompt is displayed in Figure 124.




Figure 124

Figure 125 shows the result of an EDD Data File history query. The listing show the DRG number, if the file was inserted into the database and the date as well as other general information about the file.

The screenshot shows a Netscape browser window with the title 'PEMS - Netscape'. The address bar displays the URL: https://www.pems.bechteljacobs.org/qpems/qpems/Pems.load_menus. The page content is titled 'History of Data File: R0023243.RES' and includes a red button labeled 'CLICK HERE TO CONTINUE'. A table titled 'Current File Information' displays the following data:

DRG No.	Inserted	Date Inserted	If Verified	General Info
L0016624	Inserted	03/21/2003 05:47 pm	UNLOCK	Errors: No Errors File Type: Rad Results Lab: DCSL Original Name: C:\Working\r0023243.res Project: INSTMAN Statement Of Work: BKET01320 Userid: QA
L0016623	Not loaded			Errors: Errors File Type: Rad Results Lab: THEROR Original Name: C:\Working\r0023243.res Project: INSTMAN Userid: QA EDD Date: 03/21/2003 05:45 pm
L0016622	Not loaded			Errors: Errors File Type: Rad Results Lab: THEROR Original Name: C:\Working\r0023243.res Project: INSTMAN Userid: QA EDD Date: 03/21/2003 05:43 pm
L0016621	Not loaded			Errors: Errors File Type: Rad Results

The sidebar on the left contains the following links: Instruction Manual Project, Lab Data, Data Management, Global Admin, Project Admin, Reports, What's New, Switch Project, and Help.

Figure 125

6.4 SOW Detail Information

The SOW Detail report shows a listing of all samples by SOW and their status. It assists the Laboratory in knowing whether or not the project has finalized the LCOC so that the data can be loaded. By clicking on 'SOW Detail Information', the user can review general SOW information by SOW number and sample id. Figure 126 shows an example of this feature.

SOW Detail Information

SOW	Laboratory Code	Project ID	LCOC	LCOC Finalized?	Client Sample Id	Sample Date	Shipped Date	DRG No.
BKET01320	DCSL	INSTMAN	INST_1	Y	TSKREFADD1	03/21/2003	03/21/2003	

Figure 126

6.5 Missing Analytes Report

The 'Missing Analyte Report' (Figure 127) gives the user an overview of any missing analytes for a particular Parameter Analysis Group for a particular project. This is helpful for large projects with samples having several methods and multiple analytes being reported. This report is limited to those SOWs in Tracker.




Figure 127

As a safeguard, PEMS gives a unique database name (DRG Number) each time an attempt is made to load the file, as shown in Figure 128. Resubmission is triggered when a user loads a file name that matches an already loaded file name. For example, if 'r1234567.res' is loaded successfully and then the same name is used again, the resubmission warning is triggered. The DRG Number it is used to identify the data that must be deleted. Therefore, if new data is to be added to a data file, the new file must contain all of the previous data plus the data to be added (Figure 129).




Figure 128




Figure 129

7.0 Analytical Data

7.1 Validation Qualifier Assignment

If the analytical data is to be validated, the following applications may be used to assign validation qualifiers and reason codes that explain why the validation qualifier was used. Lab qualifiers that match existing validation qualifiers can be copied to the validation qualifier field. Sample results can be edited using the 'Edit Sample Result' button. The user can also copy and edit the qualifiers and edit sample results after resubmission. The following figures illustrate these features. Click on the Data Management menu button, then 'Validation Qualifier Assignment' (Figure 130) under "Analytical Data."

The screenshot shows a Netscape browser window titled "PEMS - Netscape". The address bar displays the URL: https://www.pems.bechteljacobs.org/qpems/qpems/Pems.load_menus. The page content is titled "Validation Qualifier Assignment". On the left, there is a vertical navigation menu with the following items: "Instruction Manual Project", "Lab Data", "Data Management" (highlighted in yellow), "Global Admin", "Project Admin", "Reports", "What's New", "Switch Project", and "Help". The main content area contains a form with the following fields and buttons:

Laboratory:	All Labs	
Filename:		Search
Sample ID:		Search
Analysis Type:		Search

Below the form, there are four buttons arranged in two rows:

- Copy Lab Qualifiers
- Edit Validation Qualifiers
- Assign A Single Qualifier
- Edit Sample Results

At the bottom, there is a single button: Assign Qualifier to Unvalidated Results.

Figure 130

Clicking on the ‘Copy Lab Qualifiers’ button as shown in Figure 131, will copy the lab qualifier assigned by the laboratory in the AMSED EDD to the validation qualifier field. Laboratory results that are not reported with a result qualifier will be assigned the “=” qualifier if the “Copy” function is used.

Only Lab Qualifiers that match OREIS validation qualifiers are subject to the copy function. Blank Lab Qualifiers will be converted to the " = " validation qualifier

Copy Sample Result Qualifier

Project: INSTMAN
 Lab Code: All File Name : All
 Sample ID: All Analysis Type : All

	Lab Qualifier	Qualifier Class	Sample ID	Analysis Method ID	Parameter	Result	Unit	Use Flag	Replicate Number	Validation Qualifier
<input checked="" type="checkbox"/>			TSKREFADD2	RAD-TC-99 BY LSC	Technetium-99	297.26	PC/L	Y	0	=
<input checked="" type="checkbox"/>			TSKREFADD2	RAD-TOTAL SR BY BETA	Radioactive Strontium (Total)	7.88	PC/L	Y	0	=
<input checked="" type="checkbox"/>			TSKREFADD2	RAD-U ISO BY ALPHA	Uranium-233/234	42.87	PC/L	Y	0	U
<input checked="" type="checkbox"/>			TSKREFADD2	RAD-U ISO BY ALPHA	Uranium-235	14.25	PC/L	Y	0	=
<input checked="" type="checkbox"/>			TSKREFADD2	RAD-U ISO BY ALPHA	Uranium-236	1.68	PC/L	Y	0	=
<input checked="" type="checkbox"/>			TSKREFADD2	RAD-U ISO BY ALPHA	Uranium-238	100.91	PC/L	Y	0	=

Copy Lab Qualifier

Figure 131

The ‘Assign Valid Qualifier and Reason Code’ screen allows the user to assign a single validation qualifier to all selected analytical results and associated reason codes if required. Valid validation qualifiers and reason codes are available in a drop down menu. Click on the appropriate qualifier and reason and then the “ok” button to apply them as shown in Figure 132. It should be noted that this function will overwrite any qualifier that may have already been assigned.




Figure 132

The 'Edit Validation Qualifiers' screen as shown in Figure 133, allows validation qualifiers and reason codes to be applied or edited to individual analyte results. This screen allows the setting of the Use Flag, which indicates whether or not the result will be transferred to OREIS.




Figure 133

A separate browser window that contains the valid validation qualifiers is launched from the 'Edit Validation Qualifiers' screen as shown in Figure 134. If the data will go to OREIS, but is not validated, the qualifier of either "XV" or "XZ" is entered before submitting the data to OREIS.




Figure 134

Figure 135 shows separate browser window that contains the list of valid reason codes that is launched from the 'Edit Validation Qualifiers' screen.

Reason Codes

Reason Code	Description
A01	Extraction holding times were exceeded.
A02	Extraction holding times were grossly exceeded.
A03	Analysis holding times were exceeded.
A04	Analysis holding times were grossly exceeded.
A05	Samples were not preserved properly.
A06	Professional judgement was used to qualify the data.
B01	Mass calibration was in error, even after applying expanded criteria.
B02	Mass calibration was not performed every 12 hours.
B03	Mass calibration did not meet ion abundance criteria.
B04	Professional judgement was used to qualify the data.
B07	Blank data not reported. Validation incomplete.
C01	Initial calibration RRF was <0.05.
C02	Initial calibration RSD was >30%.
C03	Initial calibration sequence was not followed as required.
C04	Continuing calibration RRF was <0.05.
C05	Continuing calibration %D was >25%.
C06	Continuing calibration was not performed at the required frequency.
C07	Resolution criteria were not met.
C08	RPD criteria were not met.
C09	RSD criteria were not met.
C10	Retention time of compounds was outside windows.

Figure 135

As shown in Figure 136, if an invalid reason code is entered, the user will be prompted to provide a valid code.




Figure 136

From the 'Edit Validation Qualifier' screen, laboratory qualifiers can be changed, added or deleted by clicking on the lab qualifier field for a sample result (Figure 137). The valid laboratory qualifiers for the selected analysis class (i.e., radiological, organic, wet chemical, physical, etc.), and the numeric result and result units will be shown and can be updated using this screen.




Figure 137

From the ‘Validation Qualifier Assignment’ screen and selecting “Edit Sample Results”, information may be updated and comments may be added by clicking on the record to be updated (Figure 138).




Figure 138

7.2 Edit Sample Results

The 'Edit Sample Results' screen, Figure 139, allows the user to make general and revision comments. The following fields can be updated or changed via this screen: Result value, result unit, analysis type, counting error, dilution factor, and use flag.

PEMS - Netscape

File Edit View Go Communicator Help

Back Forward Reload Home Search Netscape Print Security Shop Stop

Location: https://www.pems.bechteljacobs.org/qpems/qpems/Pems.load_menus

Instant Message Internet EIM Sites Lookup New&Cool LMES Internal H Oracle Document CNET.com - Welc WebWHOS: Search Seagull's Therm

Instruction Manual Project

Lab Data

Data Management

Global Admin

Project Admin

Reports

What's New

Switch Project

Help

Edit Customer Sample Result

Required fields are marked with *

Project ID:	INSTMAN
Lab Code:	DCSL
File Name:	R0023243.RES
Sample ID:	TSKREFADD1
Parameter:	Technetium-99
Anly Method ID:	RAD-TC-99 BY LSC
Lab Sample ID:	03-02093-04
Param Anly Group:	1227CHEM1
Replicate Number:	0
Lab Detection Limit:	4.4
Revision Flag:	N
General Comment:	<input type="text"/>
Revision Comment:	<input type="text"/>
Result Value:	<input type="text" value="297.26"/>
Result Unit:	<input type="text" value="PCI/L"/>
Analysis Type:	<input type="text" value="RADS"/>
Counting Error:	<input type="text" value="5.2"/>
Dilution Factor:	<input type="text" value="1"/>
Use Flag:	<input type="text" value="v"/>

Document: Done

Figure 139

7.3 Apply Data Assessment Codes

The Apply Data Assessment Codes function, as seen in Figure 140, applies data assessment codes and reason codes to regular analytical results and updates the use flag for these results. The function prompts the user for the name of the file containing the result record identifier, the data assessment code, the reason codes (if any) and the use flag. This file is in a specified format, which can be reviewed on the PEMS Home Page.

On clicking the Continue icon, the file is examined for invalid record identifiers, duplicated record identifiers, invalid data assessment codes, invalid reason codes and other invalid conditions.




Figure 140

If any errors are detected in the validation process, they will be displayed to the user as shown in Figure 141. This screen lists the record number in the data file that contained the error, the value of the field in error and a descriptive message of the problem encountered. These errors must be resolved by correcting the submitted file. Clicking on the 'Return to Main Screen' button will return the user to the file prompt (Figure 140).




Figure 141

If there are no errors in the file, a screen with a successful message (See Figure 142) is displayed. Clicking on the 'Apply Codes to Results' button will apply the values listed in the data assessment file to the analytical result records. If a validation qualifier is not already assigned to these records, the validation qualifier will be updated to 'DA'.




Figure 142

After the analytical result records have been updated, a screen will be displayed that reports the number of records affected by the update and a unique file identifier for the data assessment file (See Figure 143). Clicking on the 'Ok' button will return the user to Data Management menu.




Figure 143

Data Assessment codes and their associated reason codes can be viewed on the Edit Validation Qualifiers screen (see Figure 133). They can also be downloaded through the Download Sample Measurement Report.

7.4 Result/DRG Verification

The Result/DRG verification function as seen in Figure 144, allows Data Receipt Groups (DRG) from the selected laboratory to be verified or unverified. Verifying a DRG prevents resubmission of the file (Figure 145). Un-verifying a DRG will allow resubmission of the file.




Figure 144




Figure 145

The project or data validation user can get a report of validation qualifiers, QC, TIC, and TIR result information. Its purpose is to display the validation information of analytical results and their Reason Codes. This information can be displayed on the screen or downloads the data to a CSV file on the user's PC, which can be opened with Excel (Figure 146).

Validation Qualifier Report

Project ID: INSTMAN
 Lab Code: All File Name: R0023243.RES
 Sample ID: All Analysis Type: All

Sample ID	Analysis Method	Lab Sample ID	Parameter/Description	Sample Result	Lab Qualifier	Validation Qualifier	Reason Code(s)	Units	Error	Detect Limit	TPU	TPU Level
TSKREFADD1	RAD-TC-99 BY LSC	03-02093-04	14133767 Technetium-99	297.26000				PC/L	5.2	4.4	45.87	2
TSKREFADD1	RAD-TOTAL SR BY BETA	03-02093-04	NS951 Radioactive Strontium (Total)	7.88000				PC/L	1.93	1.37		
TSKREFADD1	RAD-U ISO BY ALPHA	03-02093-04	NS632 Uranium-233/234	42.87000				PC/L	8.42	.37	8.87	2
TSKREFADD1	RAD-U ISO BY ALPHA	03-02093-04	15117961 Uranium-235	4.25000				PC/L	1.83	.46	1.85	2
TSKREFADD1	RAD-U ISO BY ALPHA	03-02093-04	13982702 Uranium-236	1.68000				PC/L	1.05	.41	1.05	2
TSKREFADD1	RAD-U ISO BY ALPHA	03-02093-04	24678828 Uranium-238	67.91000				PC/L	12.49	.64	13.24	2

CLICK HERE TO CONTINUE

OR

Download to Excel

Figure 146

8.0 RTL Review and Generation

8.1 RTL Review

The RTL Review function allows the user to see any gaps in the information to be sent to OREIS (Figure 147). This report performs almost all of the validation that OREIS performs on a received data package.

The screenshot shows a Netscape browser window titled "PEMS - Netscape" with the address bar displaying "https://www.pems.bechteljacobs.org/qpems/qpems/Pems.load_menus". The page content is titled "PEMS to OREIS RTL Review Main Screen" and includes a sub-header "This module checks PEMS data for OREIS RTL errors". A note states: "Please note that the minimum and maximum possible sample collection/field measurement dates for this project are: 01/01/1992 through 04/20/2004". Below this, it says "Required fields are marked with" followed by a red circle with a white arrow icon. The form contains the following fields:

Project Task:	<input type="text" value="All Tasks"/>
Start Sample Collection Date(mm/dd/yyyy):	<input type="text"/>
End Sample Collection Date(mm/dd/yyyy):	<input type="text"/>
RTL Destination	
<input checked="" type="radio"/>	Oak Ridge OREIS
<input type="radio"/>	Portsmouth OREIS

At the bottom of the form, there is a red button labeled "CLICK HERE TO CONTINUE". The left sidebar of the browser contains a menu with the following items: Instruction Manual Project, Lab Data, Data Management (highlighted), Global Admin, Project Admin, Reports, What's New, Switch Project, and Help. The Windows taskbar at the bottom shows the Start button, several open application windows (M..., A..., P..., R..., P...), and the system clock indicating 9:16 AM.

Figure 147

Note the missing data in Figure 148. After missing data is entered into PEMS, the data is ready to send to OREIS.

PEMS - Netscape

File Edit View Go Communicator Help

Back Forward Reload Home Search Netscape Print Security Shop Stop

Bookmarks Location: https://198.124.46.174/pems/ppems/Pems.load_menus What's Related

Instant Message WebMail Calendar Radio People Yellow Pages Download Customize...

PEMS to OREIS RTL Review

SAMP MEAS Mandatory fields with no values

Field Name	Sample ID	Parameter
MED_TYPE	EGGS1	15438310
	EGGS1	N704

Field Meas Mandatory Fields With No Values

Field Name	Station	Date Collected	Parameter
MED_TYPE	0125	04/07/2003	N704
	0125	04/07/2003	NS998

CLICK HERE TO CONTINUE

Document: Done

Start Micro... Ameri... PE... PEM... 7:35 AM

Figure 148

8.2 Send a RTL to OREIS

Follow the instructions on the screen (Figure 149) and complete the 'OREIS Environmental Measurements Data Transmittal' form and return it to OREIS staff.




Figure 149

The user selects the date range to be transferred to OREIS as shown in Figure 150.

The screenshot shows a Netscape browser window titled "PEMS - Netscape". The address bar displays the URL: https://www-pems.bechteljacobs.org/qpems/qpems/Pems.load_menus. The page content includes a navigation menu on the left with links: Instruction Manual Project, Lab Data, Data Management (highlighted), Global Admin, Project Admin, Reports, What's New, Switch Project, and Help. The main content area features a header with contact information for Teresa James, Environmental Information Manager for Bechtel Jacobs, LLC. Below this, a note specifies the minimum and maximum sample collection/field measurement dates as 01/01/1992 through 04/20/2004. A section titled "Required fields are marked with" followed by a red circle with a white dot. The form contains a "Project Task" dropdown menu set to "All Tasks", two input fields for "Start Sample Collection Date(mm/dd/yyyy)" and "End Sample Collection Date(mm/dd/yyyy)", each with a red circle with a white dot indicating it is a required field. Below these is a section titled "RTL Destination" with two radio button options: "Oak Ridge OREIS" (selected) and "Portsmouth OREIS". A note at the bottom states: "Note: This transfer excludes analytical data from Tracer and Surrogate analysis." A large red button labeled "CLICK HERE TO CONTINUE" is positioned at the bottom of the form area. The Windows taskbar at the bottom shows the Start button, several open application windows, and the system clock displaying 9:20 AM.

PEMS - Netscape

File Edit View Go Communicator Help

Back Forward Reload Home Search Netscape Print Security Shop Stop

Bookmarks Location: https://www-pems.bechteljacobs.org/qpems/qpems/Pems.load_menus What's Related

Instant Message WebMail Calendar Radio People Yellow Pages Download Customize...

Instruction Manual Project

Lab Data

Data Management

Global Admin

Project Admin

Reports

What's New

Switch Project

Help

For questions concerning PEMS or OREIS policy or procedures, contact Teresa James, Environmental Information Manager for Bechtel Jacobs, LLC at 865-576-2057.

Please note that the minimum and maximum possible sample collection/field measurement dates for this project are: 01/01/1992 through 04/20/2004

Required fields are marked with

Project Task: All Tasks

Start Sample Collection Date(mm/dd/yyyy):

End Sample Collection Date(mm/dd/yyyy):

RTL Destination

☒ Oak Ridge OREIS

☐ Portsmouth OREIS

Note: This transfer excludes analytical data from Tracer and Surrogate analysis.

CLICK HERE TO CONTINUE

Document: Done

Start M... A... P... R... P... 9:20 AM

Figure 150

The user selects the date range to be transferred to OREIS. The system will ask the user to confirm the information. (Figure 151).




Figure 151

The system will rerun the review application and will check for missing data as shown in Figure 152.




Figure 152

The system will confirm that the data set has been sent to OREIS as shown in Figure 153. At this time the user may download any of the data sets that were sent to OREIS to a PC if desired.




Figure 153

This form must be sent to the OREIS support personnel listed in the instructions.

OREIS Environmental Measurements Data Transmittal

I. CONTACTS			
1. Data Custodian/Organization		2. Phone	3. Date
4. Data Generator/Organization		5. Phone	
II. PROJECT INFORMATION			
6. Project Name		7. Project Code	
8. Date Project Initiated (mm/dd/yy)		9. Date Project Completed (mm/dd/yy)	
10. ADS Name		11. ADS Number	
12. Project Sponsor <input type="radio"/> ER <input type="radio"/> EM <input type="radio"/> Other (please specify)		13. Project Program	
14. Site <input type="radio"/> ORNL <input type="radio"/> K-25 <input type="radio"/> Y-12 <input type="radio"/> Paducah <input type="radio"/> Portsmouth <input type="radio"/> Offsite <input type="radio"/> Other (specify)			
15. Project Records Location (complete location information for lab books, data packages, QA/QC, etc.)			
16. Project Description			
III. DATA INFORMATION (attach pages as needed)			
17. General Description of Data			
18. Data Files	File/Table	# Records	Metadata
19. Data Format <input type="radio"/> TXT <input type="radio"/> DBF <input type="radio"/> XLS <input type="radio"/> MDB <input type="radio"/> CSV <input type="radio"/> Other (specify)			
20. Data Transfer Media Description			
21. Footnote (brief statement that could be used with tables or graphs)			
22. Time Period (beginning yr. - ending yr.)			
23. Cautions Associated with Using the Data			
24. Complete References to Relevant Reports or Documents			
25. Exact Location of Information in Support of this Transmittal (if not included in #15)			
IV. OREIS Receiving/Review			
26. OREIS Login By (name of OREIS Data Management Staff Member)	27. Phone	28. Date Received	29. Date Acknowledged

30. Comments	31. PROJECT_ID	32. TRANS_ID	33. Document/Data Clearance Form Included

9.0 Utility Functions

9.1 Delete/Rename Sample

For a variety of reasons sample IDs may need to be edited. The Delete/Rename Sample function allows users to delete samples or to change the sample ID (rename) (Figure 154).




Figure 154

Figures 155 and 156 illustrate renaming a sample.

The screenshot shows a Netscape browser window titled "PEMS - Netscape". The address bar displays the URL "https://www.pems.bechteljacobs.org/qpems/qpems/Pems.load_menus". The browser's toolbar includes buttons for Back, Forward, Reload, Home, Search, Netscape, Print, Security, Shop, and Stop. Below the toolbar is a bookmarks bar with various links. On the left side of the browser window is a vertical navigation menu with the following items: Instruction Manual Project, Lab Data, Data Management (highlighted in yellow), Global Admin, Project Admin, Reports, What's New, Switch Project, and Help. The main content area of the browser displays a form titled "Enter the New Sample ID". The form contains three input fields: "Project ID:" with the value "INSTMAN", "Old Sample ID:" with the value "WB-01", and "New Sample ID:" which is currently empty. Below these fields are two buttons: "OK" and "Back to Query". The browser's status bar at the bottom shows "Document: Done".

PEMS - Netscape

File Edit View Go Communicator Help

Back Forward Reload Home Search Netscape Print Security Shop Stop

Bookmarks Location: https://www.pems.bechteljacobs.org/qpems/qpems/Pems.load_menus What's Related

Instant Message Internet EIM Sites Lookup New&Cool LMES Internal H Oracle Document CNET.com - Welc WebWHOS: Search Seagull's Therm

Instruction Manual Project

Lab Data

Data Management

Global Admin

Project Admin

Reports

What's New

Switch Project

Help

Enter the New Sample ID

Project ID: INSTMAN

Old Sample ID: WB-01

New Sample ID:

OK Back to Query

Document: Done

Figure 155




Figure 156

Figure 157 shows the message displayed by PEMS when a sample ID is deleted.




Figure 157

9.2 Change Sample Status

PEMS allows projects the ability to change the status of a sample. A sample can be 'retired' or marked as 'uncollectible' as shown in Figure 158.




Figure 158

9.3 Change Monitoring Location

Figure 159 shows the screen from which a monitoring location may be changed.




Figure 159

10.0 Reports

Figure 160 displays the Reports menu. A table with a short description of each report follows.




Figure 160

Report Name	Purpose
Project Report	List of projects, description and start date
PEMS Report	By project and list of fields to pick and include
Analytical Data Package Review Checklist	By LCOC from list, can print for validation
Validation Qualifier Assignment Report	By laboratory, file name, sample id and analysis type
Sample Query	By sample id for a project, gives sample information/results
Sample Status	By status code and date range
Project Status Report	By site code
Valid Value Report	Takes the user to Global Administration
Download Sample Measurements	By date range, DRG #, location and task code
Download Field Measurements	By date range, location and task code
Download Associate Samples	
Download Station	Used to check the RTL before sending data to OREIS
Download Associate QC Samples	Download to Excel reg. samples with QC
Generate DMR DBF File	By month and year
Method Blank Report	By pick list of file names
Download Sample Tracking	By date range, pick list of projects
PEMS to OREIS RTL Review	By Task Code and Sample Collection Date range, perform OREIS checks on RTL data
Download QC Data	By Analysis Date range, DRG Number, Sample Delivery Group and/or Batch Number, download QC analytical result data (LCS, Method Blank, etc.)
Sample Result and Location Report	By Sample Collection Date range, SOW, DRG Number, Monitoring Location, and/or Task Code, download a user-defined set of data fields for TIC, TIR and regular analytical results
EDD Review Reports	Purpose
QC File Confirmation Report	QC files for current or all projects
Holding Time Confirmation Report	By project, pulled from Tracker
Missing Analytes Report	By parameter analysis group
Turnaround Time Report	By SOW, pulled from Tracker